




Applied Thermal Control Ltd  
39 Hayhill Industrial Estate  
Barrow-upon-Soar, Loughborough  
LE12 8LD, United Kingdom  
+44 (0) 1530 839 998  
Service@thermalexchange.co.uk  
Support@app-therm.com

# Installation, Operation & Service Manual

## R10 / R20

### DOCUMENT DETAILS

Date	9/MAR/2022	Compiled by	MJH	Page	1 / 6	Revision	1
------	------------	-------------	-----	------	-------	----------	---





**Applied Thermal Control Ltd**  
39 Hayhill Industrial Estate  
Barrow-upon-Soar, Loughborough  
LE12 8LD, United Kingdom  
+44 (0) 1530 839 998  
Service@thermalexchange.co.uk  
Support@app-therm.com

# Installation, Operation & Service Manual

## R10 / R20

### DOCUMENT DETAILS

Date	9/MAR/2022	Compiled by	MJH	Page	2 / 6	Revision	1
------	------------	-------------	-----	------	-------	----------	---

### CHANGE LOG

Date	Revision	Page ref	Change
9/MAR/2022	1	All	First release



Applied Thermal Control Ltd  
 39 Hayhill Industrial Estate  
 Barrow-upon-Soar, Loughborough  
 LE12 8LD, United Kingdom  
 +44 (0) 1530 839 998  
 Service@thermalexchange.co.uk  
 Support@app-therm.com


# Installation, Operation & Service Manual R10 / R20

## DOCUMENT DETAILS

Date	9/MAR/2022	Compiled by	MJH	Page	3 / 6	Revision	1
------	------------	-------------	-----	------	-------	----------	---

## PRODUCT SPECIFICATIONS

Attribute	R10	R20
Weight	36kg	48kg
	79.4lbs	105.8lbs
Heat transfer fluid volume	4000mL	7500mL
	1.057USGal	1.981USGal
Cooling capacity <small>(T<sub>setpoint</sub> +20°C / +68°F)</small> <small>(T<sub>ambient</sub> +20°C / +68°F)</small> <small>*(T<sub>ambient</sub> +30°C / +86°F)</small>	*1150W	*1900W
	*3924BTU/h	*6483BTU/h
	*0.327TR	*0.54TR
Temperature stability	0.1±°C	0.1±°C
Temperature resolution	0.1±°C	0.1±°C
Settable temp. range	+4 to +35°C	+4 to +35°C
Optional temp. range	0 to +65°C	0 to +65°C
Power supply requirement	230Vac	230Vac
	50Hz	50Hz
	1~	1~
	6.5A@230Vac	7.5A@230Vac
Sound pressure level	63dBa	63dBa
Controller screen size	n/a	n/a
Dimensions (W*D*H)	6U 19" rack	8U 19" rack
	560mm / 22.05" deep	560mm / 22.05" deep
Diagnostic functions	PID refrigeration system utilization, level OK lamp, cooling output active lamp.	PID refrigeration system utilization, level OK lamp, cooling output active lamp.
Fluid fittings (standard)	1/2" BSPPF	1/2" BSPPF
Fluid fittings (option)	Pushfit 12mm	Pushfit 12mm
Low fluid level alarm	Pump interlock for empty	Pump interlock for empty
Controller	Two-row display IP65-rated	Two-row display IP65-rated
Tool-less access	No	No
Overtemperature protection	Pump motor thermal c/o	Pump motor cutout
	Optional HT trip available	
Compressor overload protection	HP switch	HP switch
Overcurrent protection	Fused	Fused
Rated duty cycle	Continuous	Continuous
Compatible heat transfer fluids	DI water, propylene glycol mixes, Hexid A4 & A6	DI water, propylene glycol mixes, Hexid A4 & A6





**Applied Thermal Control Ltd**  
39 Hayhill Industrial Estate  
Barrow-upon-Soar, Loughborough  
LE12 8LD, United Kingdom  
+44 (0) 1530 839 998  
Service@thermalexchange.co.uk  
Support@app-therm.com

# Installation, Operation & Service Manual

## R10 / R20

### DOCUMENT DETAILS

Date	9/MAR/2022	Compiled by	MJH	Page	4 / 6	Revision	1
------	------------	-------------	-----	------	-------	----------	---

### PERFORMANCE CURVES

**R10**

**R20**

**TBC**

**TBC**

**TBC**

**TBC**

**TBC**

**TBC**



### DOCUMENT DETAILS

Date	9/MAR/2022	Compiled by	MJH	Page	5 / 6	Revision	1
------	------------	-------------	-----	------	-------	----------	---

### SAFETY NOTICES

For your safety, we draw your attention to the following warning and caution marks throughout the manual; the safe operation of an ATC chiller always remains the responsibility of the operator. This equipment is intended to be used as a liquid temperature conditioning device – it requires no external pump, nor any further manipulation of temperature. If the equipment is used in a manner not specified by the manufacturer, the protection provided by the equipment may be impaired.



Caution; Failure to comply with a caution will invalidate product warranty and absolve ATC from any liability, howsoever caused, and could result in permanent damage to equipment.



Caution; Filling/topping up of the tank should only be undertaken with the unit switched off, to prevent back-filling of the fluid.



Caution; The high integrity refrigeration system contains no user-serviceable parts. Repair and service requires specialized knowledge and tools to be provided by ATC or its local agent. Any unauthorized tampering with the refrigeration system automatically invalidates warranty.



Warning; Very cold surfaces and gases, lower than -20°C (-40°F). Severe frostbite hazard.



Warning; Opening the refrigeration system may expose the operator to toxic and corrosive compounds (HFCs). Take protective measures including suitable eye protection.



Warning; Gases may exceed 15 barg (220psig) during operation.



Warning; Refrigerant is class A1. It does not support combustion but is oxygen depleting. Review equipment rating label for specific refrigerants and CO<sub>2</sub>e.



Warning; Water and electricity are in close proximity. Always ensure the unit is isolated before service. The R-series is protected from overcurrent by mains fusing. Never bypass these components.



Warning; Failure to comply with a 'warning' may result in personal injury or death. ATC does not accept any liability for injury caused through use of this equipment.



Warning; After switching off, the fan blades slow to a stop. Do not open until the fan has stopped rotating.



**Applied Thermal Control Ltd**  
 39 Hayhill Industrial Estate  
 Barrow-upon-Soar, Loughborough  
 LE12 8LD, United Kingdom  
 +44 (0) 1530 839 998  
 Service@thermalexchange.co.uk  
 Support@app-therm.com

# Installation, Operation & Service Manual

## R10 / R20

### DOCUMENT DETAILS

Date	9/MAR/2022	Compiled by	MJH	Page	6 / 6	Revision	1
------	------------	-------------	-----	------	-------	----------	---

### INCLUDED ANNEXES

Specific technical product information is provided in the following series of annexes.

- Annex A-1 191122 Shipping & Unpacking - Weighing over 18kg without castors
- Annex B-13 220309 Site & Environmental Requirements for R10 & R20
- Annex C-2 220309 Installation - Generic water-cooled with 0.5inch fittings
- Annex C-6 220309 Installation - Generic air-cooled with 0.5inch fittings
- Annex D-3 201007 Fluid Handling & Startup Procedures - Generic
- Annex E-15 210301 KR1 - how to use
- Annex E-15C 210916 KR1 Program General Use K1 thru K12
- Annex F-3 201007 PD pump using discrete PRV
- Annex G-8 220309 Troubleshooting - Initial help air-cooled K-, R- & G-Series
- Annex G-14 220309 Troubleshooting - Initial help water-cooled K-, R- & G-Series
- Annex H-1 191121 End-user maintenance - air-cooled units with water as fluid
- Annex H-2 210326 End-user maintenance - water-cooled units with water as fluid
- Annex I-1 210830 Maintenance for technicians - Generic refrigerated units
- Annex J-5 200706 EU Compliance Statement - Conflict Minerals
- Annex J-7 200715 EU Compliance Statement - REACH
- Annex J-8 200827 EU Compliance Statement - POPs
- Annex J-10 201111 EU Compliance Statement - RoHS
- Annex J-12 210830 UKCA DoC - R-series
- Annex J-13 210830 EU DoC - R-series
- Annex K-1 200623 Standard warranty terms of ATC
- Annex M-9 220309 Recommended spares, R10 & R20
- Annex R-1 170621 SDS Refrigerant HFC-R134a
- Annex R-2 120701 SDS Refrigerant HFC-R407C
- Annex R-3 200203 SDS Hexid A4 v6.4



Applied Thermal Control Ltd  
39 Hayhill Industrial Estate  
Barrow-upon-Soar, Loughborough  
LE12 8LD, United Kingdom  
+44 (0) 1530 839 998  
Service@thermalexchange.co.uk  
Support@app-therm.com

# Annex A-1

## DOCUMENT DETAILS

Date	22/NOV/19	Author(s)	MJH	Page	1 / 1	Revision	2
------	-----------	-----------	-----	------	-------	----------	---

### UNPACKING UNITS WEIGHING OVER 18kg (40lbs) WITHOUT CASTORS

1	Please check that both the packaging and the unit are undamaged. If there is any doubt, it is vital that you inform both ATC and the carrier. There are no hidden shipping bolts or other fixings. You should inspect the packaging for signs of transit damage before signing for the unit, and if possible, unpack the unit before signing. Once you have signed for the goods, ATC cannot be held responsible for any transit damage subsequently found.
2	As the unit is >18kg, ATC must recommend that 2 persons are used to lift by hand, or a crane.
3	Remove the unit from its original packaging and ensure that there is no packaging left around the cooling ducts. There is no internal product packaging that requires the chiller to be opened.
4	<b>Please retain all packaging in the unlikely event that the chiller needs to be returned to our local representatives.</b>



Applied Thermal Control Ltd  
 39 Hayhill Industrial Estate  
 Barrow-upon-Soar, Loughborough  
 LE12 8LD, United Kingdom  
 +44 (0) 1530 839 998  
 Service@thermalexchange.co.uk  
 Support@app-therm.com

# Annex B-13

## DOCUMENT DETAILS

Doc. Date	9/MAR/2022	Author(s)	MJH	Page	1 / 1	Revision	1
-----------	------------	-----------	-----	------	-------	----------	---

## R10 & R20 SITE & ENVIRONMENTAL REQUIREMENTS

**1 Rack mount** – this product is only intended to be used with a ventilated standard 19” rack cabinet. R10 is a 6U high unit and R20 is an 8U high unit. Both units are 560mm deep. Both units require 4-off end user-supplied M6 bolts to secure to rack cabinet upright rails. ATC strongly recommends the use of support guides that run from the front upright rail to the back upright rail of the rack cabinet to support the weight of the chiller.

**2 Clean, dust-free environment** – the R10 & R20 units use external axial frame fans to provide airflow for cooling the pump motor, compressor and condenser. Do not block these vents. The system will shut down in the event the motor high temperature or refrigeration system high pressure switches are triggered.

**3 Non-condensing ambient temperature** – +5°C to +40°C (+39°F to +104°F). This prevents build-up of moisture on internal components.

**4 Humidity** - 80% for ambient temperatures up to +31°C (+88°F), decreasing linearly to 50% relative humidity at +40°C (+104°F) ambient temperature.

	Electrical supply	R10		R20	
	ATC power supply specification	-0spec	-9spec	-0spec	-9spec
Nominal supply voltage	230Vac	208-230Vac	230Vac	208-230Vac	
Voltage fluctuations	±10%	±10%	±10%	±10%	
Frequency	50Hz	50/60Hz	50Hz	50/60Hz	
Mode of supply	L/N/E	L/N/E or L1/L2/E	L/N/E	L/N/E or L1/L2/E	
Current draw at 230Vac 50Hz	5A	TBC	8A	8.5A	


The power inlet module itself is rated 10A 250V. The user must provide protective earth at IEC type C14 appliance inlet. See product rating label for fuse specification. Mains fuses are typically UL-compliant and T speed rated with high breaking (rupture) capacity. Supply fuses are rear mounted. Secondary fuses protecting individual components are internally-mounted.

**6 Clearance** – the rack form factor limits clearance considerations to the rear face of the unit. Ensure there is sufficient room at the rear of the product to remove the power cable without obstruction. Ensure bend radius of your chosen hose/tube/pipework is not compromised. The front face of the R10 has no vents, but the R20 does – it is a suction vent. Left and right-hand side clearance is dictated by the rack cabinet, but both R10 & R20 have suction vents on those faces.

**7 Plumbing** – tubing, piping or hose must be clean and compatible with the fluid to be used. The chiller is compatible with deionized water and water-glycol mixtures such as Hexid fluid.

**8 Indoor use only** – altitude up to 2000m.

**9 Installation category** – transient overvoltage category II; Pollution degree 2. Temporary overvoltages occurring on mains supply are acceptable within limits defined in the aforementioned categories.

 Caution; Always use ATC recommended fluids in your chiller – many other anti-freeze mixtures have the potential to corrode your application and to damage seals in the chiller.





**DOCUMENT DETAILS**

Date	9/MAR/2022	Author(s)	MJH	Page	1 / 2	Revision	3
------	------------	-----------	-----	------	-------	----------	---

**INSTALLATION FOR WATER-COOLED UNITS WITH 1/2" BSPPF FITTINGS**

This guide applies to the following product groups;

- K-Series, R-Series and G-Series refrigerated units, where heatload is carried away by water supply.
- XR- and XF-Series water-to-water heat exchangers, where heatload is carried away by water supply.

**HOSE RECOMMENDATIONS**

Having ensured that your installation meets all site requirements, it is best practice that the fluid lines between your application and the chiller have the following characteristics.

- |   |   |
|---|---|
| 1 | <b>Short in length</b> – this reduces friction-based pressure drop and addition ambient heat load.  |
| 2 | <b>Large diameter bore</b> – at least 12mm (1/2”).  |
| 3 | <b>Free from 90° bends</b> – to limit the effects of water hammer. If this cannot be avoided, sharp changes of direction should be minimized so far as possible. Doing this correctly can yield higher pump performance and extend time between maintenance intervals. It will also reduce electrical energy consumption.   |
| 4 | <b>Clean</b> – If your installation is to existing pipe work, it is good practice to flush the system with either a commercially available central heating cleaner or 5% acetic acid solution. The system should be flushed clean with tap water to remove all traces of cleaner prior to filling the system. Failing this, it is recommended to use a domestic bleach in solution with tap water, diluted to the point where the bleach can longer be smelled by human nose. |
| 5 | <b>Opaque, ideally black</b> – to inhibit light passing through the tube and algae building up. Alternatively, solid ABS or copper pipe can be used where application chemistry allows.   |
| 6 | <b>Insulation, where low temperature process is planned</b> – the process line from chiller to application contains the feed of low temperature fluid. Insulation prevents heat from entering this line and can promote better stability. Uninsulated return lines are helpful where free cooling can be obtained by allowing heat to transfer to air – likewise, insulating the return line is helpful if the fluid temperature is below ambient.                            |



Caution; Never use transparent tubing. UV light will pass through, prompting growth of organic contamination.

**CONNECTING ADAPTERS TO PRODUCT BULKHEAD FITTINGS**

- |   |  |
|---|--|
| 1 | Standard units are supplied with 1/2” British Standard Pipe Parallel Female (BSPPF) threads (also known as G threads (ISO228)) by default, for both application/process and house/primary water supply. These fittings are not valved and will ‘drop’ the volume of the system if left open to atmosphere.         |
| 2 | Ensure the appropriate thread sealants are used in the fitting of adapters to hose. For metallic mating parts, we recommend Loctite 577. For plastic adaptors such as those supplied with the product, we recommend using ~8-12mm wide PTFE tape, wrapped around the male thread before tightening.                |
| 3 | Ensure that the system is correctly connected. The ‘donut’ labels around the ports are clearly marked with inlet and outlet symbols and function in both English and French language. Ports marked as outlets mean fluid leaves the product and must be connected to the process inlet or house water return line. |
| 4 | Check all joints are tight and leak free.  |
| 5 | Where this product is incorporated into other equipment, it is the responsibility of the assembler to ensure safety.   |



Applied Thermal Control Ltd  
39 Hayhill Industrial Estate  
Barrow-upon-Soar, Loughborough  
LE12 8LD, United Kingdom  
+44 (0) 1530 839 998  
Service@thermalexchange.co.uk  
Support@app-therm.com

# Operating Manual; Installation Annex C-2

## DOCUMENT DETAILS

Date	9/MAR/2022	Author(s)	MJH	Page	2 / 2	Revision	3
------	------------	-----------	-----	------	-------	----------	---

## BACKFILLING

- 1 In situations where the chiller is situated physically lower than the application being cooled, fluid will apply pressure to the water circuit of the product.
- 2 The weakest seal is normally the tank lid, and this is typically where fluid will escape the unit.
- 3 Ideally, the product will be located higher or level with the product water-line. If this is not possible, a non-return solenoid valve kit can be installed as an optional standard assembly.
- 4 Please raise any questions with the sales team on sales@app-therm.com.



**DOCUMENT DETAILS**

Date	9/MAR/2022	Author(s)	MJH	Page	1 / 2	Revision	2
------	------------	-----------	-----	------	-------	----------	---

**INSTALLATION FOR AIR-COOLED UNITS WITH 1/2" BSPPF FITTINGS**

This guide applies to the following product groups;

- K-Series, R-Series and G-Series refrigerated units, where heatload is rejected to air.
- A-Series airblast units, where heatload is rejected to air.

**HOSE RECOMMENDATIONS**

Having ensured that your installation meets all site requirements, it is best practice that the fluid lines between your application and the chiller have the following characteristics.

- |          |   |
|----------|---|
| <b>1</b> | <b>Short in length</b> – this reduces friction-based pressure drop and addition ambient heat load.  |
| <b>2</b> | <b>Large diameter bore</b> – at least 12mm (1/2”).  |
| <b>3</b> | <b>Free from 90° bends</b> – to limit the effects of water hammer. If this cannot be avoided, sharp changes of direction should be minimized so far as possible. Doing this correctly can yield higher pump performance and extend time between maintenance intervals. It will also reduce electrical energy consumption.   |
| <b>4</b> | <b>Clean</b> – If your installation is to existing pipe work, it is good practice to flush the system with either a commercially available central heating cleaner or 5% acetic acid solution. The system should be flushed clean with tap water to remove all traces of cleaner prior to filling the system. Failing this, it is recommended to use a domestic bleach in solution with tap water, diluted to the point where the bleach can longer be smelled by human nose. |
| <b>5</b> | <b>Opaque, ideally black</b> – to inhibit light passing through the tube and algae building up. Alternatively, solid ABS or copper pipe can be used where application chemistry allows.   |
| <b>6</b> | <b>Insulation, where low temperature process is planned</b> – the process line from chiller to application contains the feed of low temperature fluid. Insulation prevents heat from entering this line and can promote better stability. Uninsulated return lines are helpful where free cooling can be obtained by allowing heat to transfer to air – likewise, insulating the return line is helpful if the fluid temperature is below ambient.                            |



Caution; Never use transparent tubing. UV light will pass through, prompting growth of organic contamination.

**CONNECTING ADAPTERS TO PRODUCT BULKHEAD FITTINGS**

- |          |   |
|----------|---|
| <b>1</b> | Standard units are supplied with 1/2” British Standard Pipe Parallel Female (BSPPF) threads (also known as G threads (ISO228)) by default. These fittings are not valved and will ‘drop’ the volume of the system if left open to atmosphere.   |
| <b>2</b> | Ensure the appropriate thread sealants are used in the fitting of adapters to hose. For metallic mating parts, we recommend Loctite 577. For plastic adaptors such as those supplied with the product, we recommend using ~8-12mm wide PTFE tape, wrapped around the male thread before tightening. |
| <b>3</b> | Ensure that the system is correctly connected. The ‘donut’ labels around the ports are clearly marked with inlet and outlet symbols and function in both English and French language. Ports marked as outlet means fluid leaves the product and must be connected to the process inlet.             |
| <b>4</b> | Check all joints are tight and leak free.   |
| <b>5</b> | Where this product is incorporated into other equipment, it is the responsibility of the assembler to ensure safety.  |



Applied Thermal Control Ltd  
39 Hayhill Industrial Estate  
Barrow-upon-Soar, Loughborough  
LE12 8LD, United Kingdom  
+44 (0) 1530 839 998  
Service@thermalexchange.co.uk  
Support@app-therm.com

## Operating Manual; Installation

# Annex C-6

### DOCUMENT DETAILS

Date	9/MAR/2022	Author(s)	MJH	Page	2 / 2	Revision	2
------	------------	-----------	-----	------	-------	----------	---

### BACKFILLING

- 1 In situations where the chiller is situated physically lower than the application being cooled, fluid will apply pressure to the water circuit of the product.
- 2 The weakest seal is normally the tank lid, and this is typically where fluid will escape the unit.
- 3 Ideally, the product will be located higher or level with the product water-line. If this is not possible, a non-return solenoid valve kit can be installed as an optional standard assembly.
- 4 Please raise any questions with the sales team on sales@app-therm.com.



# Annex D-3

## DOCUMENT DETAILS

Date	7/OCT/2020	Author(s)	LW, MJH	Page	1 / 1	Revision	1
------	------------	-----------	---------	------	-------	----------	---

## FILLING A SYSTEM

- a) **Check all application valves are open, including solenoid valves and variable position valves.**
  - i) The product will require an open water circuit to pump into.
  - ii) Any obstructions can increase the time, or entirely prevent the bleeding of air from the system.
- b) **Remove the cap from the tank lid on the top of the product.**
  - i) Fill the tank to just underneath the bottom of the filling port.
- c) **Turn your attention to the main power switch.**
  - i) Immediately after toggling this switch assuming that the power cord is connected and turned on at the wall the product will begin to pump water.
  - ii) Leave the product to run until it cuts out on the level switch interlock.
  - iii) Turn the unit off at the mains switch.
  - iv) Fill the tank again to lift the level switches.
  - v) Start the product again until the unit cuts out again.



Caution; Do not run the pump dry. Do not deadhead the pump.

- d) **Repeat steps at c) until the chiller tank water level (if visible) doesn't drop, and the chiller doesn't cut-out on its level switches.**
- e) **With the unit now running;**
  - i) Leave the cap off the tank for >30mins to allow air to escape, or very loosely screwed on to prevent water splashing out of the tank if the unit has a 'flow through' design.
- f) **Check the application and tubing for signs of leaks whilst the chiller is running.**
  - i) Replace the tank lid fully when satisfied the system is full and bled of air.

## DRAINING A SYSTEM

- a) **Isolate the unit. Have a suitable bucket or drain on-hand.**
  - i) Remove the tank lid to allow air into the system.
  - ii) Disconnect hoses individually. Be aware hoses as well as the unit are filled with fluid.
  - iii) Consider using the red transport plugs to block product ports to give yourself time to empty hoses before continuing to empty the product.
- b) **Local rules affect where fluid can be disposed of.**
  - i) Ensure hazardous products do not enter the water course and are reclaimed from the unit for professional disposal.



# Annex E-15

## DOCUMENT DETAILS

Date	1/MAR/2021	Author(s)	MJH	Page	1 / 1	Revision	1
------	------------	-----------	-----	------	-------	----------	---

## BASIC PROGRAMMING GUIDE KR1 (ALL MODELS)

This guide may apply to your product if you require general navigation advice and help accessing settings. If you are planning to change the settings in any way, you may need a copy of the existing settings which are model dependent, signified by a letter on the end of Annex E-15, i.e. Annex E-15A.

## DISPLAY CONTENTS DURING NORMAL OPERATION

- 1 Physical navigation buttons, up, down, return and enter.
- 2 8888.8 is the actual read value on input sensor.
- 3 888.8 is the setpoint value.
- 4 Rectangles bottom left 1-4 display when output is active.
- 5 MAN LED shows in manual mode (fixed output value).
- 6 °C or °F shows units as settable in the 'inP' group.
- 7 AL LED appears when output is beyond a set alarm point.



## ACCESS TO SETTINGS

- 8 Push the return button for more than 5 seconds. The upper display will show PASS while the lower display will show 0.
- 9 Using up and down buttons set the programmed password – full access is granted by entering '40'. ATC are not responsible for damage either to the chiller or the connected equipment as a result of changing parameters without ATC's oversight.
- 10 During parameter modification the instrument continues to perform process control. In certain conditions, when a configuration change can produce a significant change to the process, it is advisable to temporarily stop the controller from controlling during the programming procedure (control outputs will be OFF). A password equal to 2000 + the programmed value (i.e. 2000 + 40 = 2040). The control will restart automatically when the configuration procedure will be manually closed.
- 11 Push the return button. If the password is correct the display will show the acronym of the first parameter group 'inP'. Push button for more than 5 seconds, the instrument will come back to the "standard display".
- 12 The configuration parameters are collected in various groups. Every group defines all parameters related with a specific function (control, alarms, output functions).
- 13 Push return button for more than 5 seconds, the instrument will come back to the "standard display". For specific settings and guidance, review the controller datasheets provided by ATC. If you are not in receipt of these, please contact ATC using the information in the header of this document.

## GENERAL NAVIGATION

- 14 Return button; A short press allows to exit from the current parameter group and select a new parameter group. A long press allows you to close the configuration parameter procedure (the instrument will come back to the "standard display").
- 15 Enter button; When the upper display is showing a group and the lower display is blank, this key allows to enter in the selected group. When the upper display is showing a parameter and the lower display is showing its value, this key allows to store the selected value for the current parameter and access the next parameter within the same group.
- 16 Up button; Allows to increase the value of the selected parameter.
- 17 Down button; Allows to decrease the value of the selected parameter.
- 18 Pushing both Return and Enter buttons moves back to the previous group. Press return first to start. The selection of the group is cyclic (on a carousel), so it is possible to move back around to the group you require.



Applied Thermal Control Ltd  
 39 Hayhill Industrial Estate  
 Barrow-upon-Soar, Loughborough  
 LE12 8LD, United Kingdom  
 +44 (0) 1530 839 998  
 Service@thermalexchange.co.uk  
 Support@app-therm.com

Operating Manual; Controller Operation

# Annex E-15C

**DOCUMENT DETAILS**

Date	16/SEP/2021	Author(s)	WE	Page	1 / 2	Revision	1
------	-------------	-----------	----	------	-------	----------	---

**STANDARD KT THRU K12 PROGRAM PURPOSE**

<b>Internal part number</b>	EA894
<b>Manufacturer part number</b>	KR1
<b>Program purpose</b>	1) STANDARD PROGRAM FOR KT THRU K12

**PROGRAM**

Group inP		Group 'out'		Group 'AL1'		Group 'AL2'	
Setting	Value	Setting	Value	Setting	Value	Setting	Value
SEnS	Pt1	o1.F	c.rEG	AL1t	LHdo	AL2t	nonE
dP	1	o1.AL	n/a	Ab1	0	Ab2	n/a
SSc	n/a	o1.Ac	dir	AL1L	-10.0	AL2L	n/a
FSc	n/a	o2F	AL	AL1H	10.0	AL2H	n/a
Unit	°c	o2.AL	1	AL1	10.0	AL2	n/a
FiL	10	o2Ac	rEU	HAL1	1.0	HAL2	n/a
inE	our	o4F	nonE	AL1d	oFF	AL2d	n/a
oPE	100	o4.AL	n/a	AL1o	0	AL2o	n/a
io4.F	out4	o4Ac	n/a				
diF1	oFF						
diF2	n/a						
di.A	0						
Group 'AL3'		Group 'LbA'		Group 'rEG'		Group 'SP'	
Setting	Value	Setting	Value	Setting	Value	Setting	Value
AL3t	nonE	LbAt	oFF	cont	Pid	nSP	n/a
Ab3	n/a	LbSt	n/a	Auto	0	SPLL	4.0
AL3L	n/a	LbAS	n/a	tunE	oFF	SPHL	35.0
AL3H	n/a	LbcA	n/a	Aut.r	n/a	SP	20
AL3	n/a			SELF	n/a	SP2	n/a
HAL3	n/a			HSEt	n/a	SP3	n/a
AL3d	n/a			cPdt	n/a	SP4	n/a
AL3o	n/a			Pb	15	A.SP	SP1
				ti	60	SP.rt	trin
				td	1	SPLr	Loc
				Fuoc	0.42	SP.u	inF
				tcH	n/a	SP.d	inF
				rcG	1.00		
				tcc	7		
				rS	0		
				Str.t	n/a		
				db.S	n/a		
				od	n/a		
				St.P	0		
				SSt	oFF		
				SS.tH	999.9		



Applied Thermal Control Ltd  
 39 Hayhill Industrial Estate  
 Barrow-upon-Soar, Loughborough  
 LE12 8LD, United Kingdom  
 +44 (0) 1530 839 998  
 Service@thermalexchange.co.uk  
 Support@app-therm.com

Operating Manual; Controller Operation

# Annex E-15C

**DOCUMENT DETAILS**

Date	16/SEP/2021	Author(s)	WE	Page	2 / 2	Revision	1
------	-------------	-----------	----	------	-------	----------	---

Group 'tin'		Group 'PrG'		Group 'PAn'		Group 'Ser'	
Setting	Value	Setting	Value	Setting	Value	Setting	Value
tr.F	nonE	Pr.F	nonE	PAS2	20	Add	1
tr.u	n/a	Pr.u	n/a	PAS3	40	bAud	9600
tr.t1	n/a	Pr.E	n/a	PAS4	300	trSP	nonE
tr.t2	n/a	Pr.Et	n/a	uSrb	tunE		
tr.St	n/a	Pr.S1	n/a	diSP	SPF		
		Pr.G1	n/a	di.CL	0		
		Pr.t1	n/a	AdE	5		
		Pr.b1	n/a	diS.t	oFF		
		Pr.E1	n/a	FiLd	0.5		
		Pr.S2	n/a	dSPu	AS.Pr		
		Pr.G2	n/a	oPr.E	ALL		
		Pr.t2	n/a	oPEr	Auto		
		Pr.b2	n/a				
		Pr.E2	n/a				
		Pr.S3	n/a				
		Pr.G3	n/a				
		Pr.t3	n/a				
		Pr.b3	n/a				
		Pr.E3	n/a				
		Pr.S4	n/a				
		Pr.G4	n/a				
		Pr.t4	n/a				
		Pr.b4	n/a				
		Pr.E4	n/a				
		Pr.St	n/a				
Group 'COn'		Group 'CAL'					
Setting	Value	Setting	Value				
Co.tY	oFF	AL.P	0				
UoLt	n/a	AL.o	!TESTRIG!				
cur	n/a	AH.P	999.9				
h.Job	n/a	AH.o	!TESTRIG!				
t.Job	n/a						




### DOCUMENT DETAILS

Date	7/OCT/2020	Author(s)	MJH	Page	1 / 2	Revision	1
------	------------	-----------	-----	------	-------	----------	---

### POSITIVE DISPLACEMENT PUMPS & DISCRETE PRESSURE RELIEF VALVE

This arrangement comprises a positive displacement pump (most commonly a rotary vane type) with a spring-loaded pressure relief valve to provide better overpressure setting control with minimal flow losses compared to fixed orifice bypasses/reliefs. This annex describes ATC's default settings and how to adjust the system.



It is important to understand the basic principle that all else being equal, higher flow results in a higher demand for pressure. The motor generates the power required to turn the pump head and create that pressure. The more restrictive a water circuit is, the higher the pressure required to maintain flowrate. Positive displacement pumps are designed to generate high pressure and are mechanically tight - their RPM dictates flowrate. See image;

- 1) **Pump discharge** – fluid leaves the pump head and enters the gauge-tee assembly. It passes the temperature sensor at the tee, and when the pressure relief is inactive, fluid heads towards no. 2.
- 2) **Outlet to process** – fluid heads out of the product to the application. Whatever restrictions lie downstream, fluid leaving this point has not passed through the pressure relief valve.
- 3) **Pressure gauge connection** – connection to a pressure gauge on the front panel of the product. This displays the pressure in the water circuit, *at the pump outlet, not the application which will be lower.*
- 4) **Pressure Relief Valve (PRV)** – spring tension in the body of the PRV determines whether fluid travels through point 2 or passes through to point 7.
- 5) **PRV adjustment knob** – rotate clockwise to increase maximum delivery pressure. Rotate anti-clockwise to limit maximum delivery pressure.
- 6) **PRV adjustment locknut** – release this locknut to adjust no.5. Ensure it is tightened once adjustments are made. If it is not, vibration may cause the adjustment knob to move on itself.
- 7) **Bypass flow outlet** – where the pressure requirement to overcome restrictions downstream of no.2 rises to be higher than the setting at no.5, the spring inside no.4 will compress and allow liquid to start bleeding through to no.7. It is important to understand that the nature of the spring means there can be no black and white point for pressure relief setting – the spring will slowly compress and bleed flow until all flow passes through no.7. When fully bypassing, all flow stays inside the chiller to protect the application.



Applied Thermal Control Ltd  
39 Hayhill Industrial Estate  
Barrow-upon-Soar, Loughborough  
LE12 8LD, United Kingdom  
+44 (0) 1530 839 998  
Service@thermalexchange.co.uk  
Support@app-therm.com

# Annex F-3

## DOCUMENT DETAILS

Date	7/OCT/2020	Author(s)	MJH	Page	2 / 2	Revision	1
------	------------	-----------	-----	------	-------	----------	---

### SETTING THE PRESSURE RELIEF VALVE (BYPASS)

Unless otherwise agreed at point of sale, the default setting for this type of pump and PRV arrangement is 3.33bar (50psi). In the absence of external calibrated pressure gauges, it's possible to use the chiller itself to set this value;

- 8) **Isolate the chiller** – always isolate before performing work.
- 9) **Disconnect from the application if already connected** – review the draining procedure in Annex D.
- 10) **Connect the chiller's process inlet to process outlet** – a short run of hose around 1-2m (3-6ft) will be adequate.
- 11) **Start the chiller and follow the fill process from Annex D** – fluid will now be running through a short loop with very low pressure required to overcome the restriction.
- 12) **With the chiller running, very slowly kink the short run of hose fitted above** – this simulates a blocked application – this is known as 'deadheading' the pump. Without PRV, the pump motor would stall or hoses might blow off. Watch the pressure gauge climb as you apply the kink.
- 13) **The pressure gauge will eventually stop climbing with the hose fully kinked** – note the value; this is the setting of the PRV where all flow is bypassed.
- 14) **Refer to points 5) and 6) above to adjust if required** – contact ATC if you're unsure over whether your desired setting is reasonable for the product you have.



# Annex G-8

## DOCUMENT DETAILS

Date	9/MAR/2022	Author(s)	MJH	Page	1 / 1	Revision	2
------	------------	-----------	-----	------	-------	----------	---

## AIR-COOLED K-, G-, R-SERIES INITIAL TROUBLESHOOTING

#	SYMPTOM	CAUSES
1	Compressor or pump motor not running Large condenser fan is running, and control illuminated.	Level switch may have been tripped – check level and top-up to resolve.
		High pressure switch may have been tripped; <ol style="list-style-type: none"> <li>Check fan spins freely and that condenser is clear of debris.</li> <li>Assess whether the ambient temperature is too high.</li> <li>Assess whether ventilation/air circulation is poor.</li> <li>Assess whether surrounding equipment is generating hot air and feeding it into the chiller air intakes.</li> </ol>
		If fitted, a low-pressure switch may have been tripped – this indicates refrigerant loss. A multimeter is required to test this.
		If running at a high setpoint, or if operating conditions have changed and raised the fluid temperature, it is possible the compressor has turned itself on high temperature internal protection.
2	Noisy operation	Pump motor or fan motor bearings may have failed.
		Pump head vane or seal may have failed.
		In a 3-phase system, phase rotation may be incorrect.
3	High fluid pressure / low flowrate	Fluid lines may have become fouled through solid debris or biological growth. Ensure biocides or pre-mixed biocidal process fluids are being used. Ensure hoses are not transparent or translucent – UV light entering prompts growth.
		General fluid leak can lead to flow rate failure.
		Excess fittings, too many sharp direction changes, too small a diameter hose for flowrate or hoses too long at the wrong diameter.
4	Fluid collecting or leaking	Leaks are impossible to diagnose remotely, but usually you will see a pool of fluid on a lower surface. Do not confuse this with condensation from cold parts of the fridge.
		Review the chemical compatibility of your fluid with the chiller. Contact ATC if you are unsure.
5	Poor cooling performance	Check fan spins freely and that condenser is clear of debris. Excess thermal load applied. Assess whether the ambient temperature is too high. Assess whether ventilation/air circulation is poor. Assess whether surrounding equipment is generating hot air and feeding it into the chiller air intakes.



# Annex G-14

## DOCUMENT DETAILS

Date	9/MAR/2022	Author(s)	MJH	Page	1 / 1	Revision	2
------	------------	-----------	-----	------	-------	----------	---

## WATER-COOLED K-, G-, R-SERIES INITIAL TROUBLESHOOTING

#	SYMPTOM	CAUSES
1	Compressor or pump motor not running whilst control remains illuminated.	Level switch may have been tripped – check level and top-up to resolve.
		High pressure switch may have been tripped; <ul style="list-style-type: none"> <li>a) Check for sufficient flowrate and pressure of house water.</li> <li>b) Assess whether the house water temperature is too high.</li> <li>c) Assess whether any fouling of condenser plate heat exchanger has occurred.</li> <li>d) Assess whether new equipment has been added to a shared water supply and contributed to either a lack of water flow (through flow splitting) or whether additional thermal load has been added to the central water's chiller/cooler.</li> </ul>
		If fitted, a low-pressure switch may have been tripped – this indicates refrigerant loss. A multimeter is required to test this.
		If running at a high setpoint, or if operating conditions have changed and raised the fluid temperature, it is possible the compressor has turned itself on high temperature internal protection.
2	Noisy operation	Pump motor bearings may have failed.
		Pump head vane or seal may have failed.
		In a 3-phase system, phase rotation may be incorrect.
3	High fluid pressure / low flowrate	Fluid lines may have become fouled through solid debris or biological growth. Ensure biocides or pre-mixed biocidal process fluids are being used. Ensure hoses are not transparent or translucent – UV light entering prompts growth.
		General fluid leak can lead to flow rate failure.
		Excess fittings, too many sharp direction changes, too small a diameter hose for flowrate or hoses too long at the wrong diameter.
4	Fluid collecting or leaking	Leaks are impossible to diagnose remotely, but usually you will see a pool of fluid on a lower surface. Do not confuse this with condensation from cold parts of the fridge.
		Review the chemical compatibility of your fluid with the chiller. Contact ATC if you are unsure.
5	Poor cooling performance	Excess thermal load applied. Assess whether the ambient temperature is too high. Assess whether house water temperature is too high. Assess whether house water flowrate is too low. Assess whether surrounding equipment is adding load or taking flow from the house water inlet to the ATC product.



Applied Thermal Control Ltd  
 39 Hayhill Industrial Estate  
 Barrow-upon-Soar, Loughborough  
 LE12 8LD, United Kingdom  
 +44 (0) 1530 839 998  
 Service@thermalexchange.co.uk  
 Support@app-therm.com

Operating Manual; Maintenance for End-Users

# Annex H-1

**DOCUMENT DETAILS**

Date November 2019

Compiled by MJH

Revision 1

**PERIODIC MAINTENANCE REQUIREMENTS BY END USER**



Caution; Failure to carry out service at the specified intervals may permanently damage your equipment.

Print this sheet out and display close to the chiller to maximize the visibility of maintenance requirements.

Weekly	Week 1	Week 2	Week 3	Week 4
Check fluid level – top up as required.				

Monthly	J	F	M	A	M	J	J	A	S	O	N	D
Check the condenser is free from dust or accumulation of debris.												

Annually	Y1	Y2	Y3	Y4	Y5	Y6	Y7	Y8
Drain process fluid and replace with fresh fluid.								
Check for fluid leaks throughout chiller and application.								
Clear any debris from inside the chiller.								

A vacuum cleaner is recommended for cleaning out the condenser, while soft cloths and IPA are recommended for cleaning metallic surfaces. If any spillages have occurred, best practice is to allow the water to evaporate off and wipe up remaining glycol residue with a cloth. Always clean with power supply isolated.



Caution; Never blow out the condenser with compressed air.



Caution; If the mains lead is lost or damaged, contact ATC who will be able to supply a replacement of the correct specification.



Applied Thermal Control Ltd  
 39 Hayhill Industrial Estate  
 Barrow-upon-Soar, Loughborough  
 LE12 8LD, United Kingdom  
 +44 (0) 1530 839 998  
 Service@thermalexchange.co.uk  
 Support@app-therm.com

# Annex H-2

## DOCUMENT DETAILS

Date	26/MAR/2021	Author(s)	MJH	Page	1 / 1	Revision	3
------	-------------	-----------	-----	------	-------	----------	---

## PERIODIC MAINTANENCE FOR WATER-COOLED PRODUCTS

This guide may apply to your unit if you have a water-condensed refrigerated unit, or a water-cooled water-to-water heat exchanger. These units do not rely on large/multiple cooling fans to reject heat to air, but instead reject heat into a house/primary/building/city water supply.

## SCHEDULE & RECOMMENDED ACTIONS

	Caution: Failure to carry out service at the specified intervals may permanently damage your equipment.
	Caution: If the mains wiring becomes damaged, contact ATC or a qualified electrician who will be able to supply a replacement of the correct specification.
	Soft cloths and IPA are recommended for cleaning metallic surfaces. If any spillages have occurred, best practice is to allow the water to evaporate off and wipe up remaining glycol residue with a cloth. Always clean with power supply isolated.
	Print this sheet out and display close to the product to maximize the visibility of maintenance requirements.

<b>Weekly</b>	<b>Check fluid level – top up as required.</b>											
	1	2	3	4	5	6	7	8	9	10	11	12
	15	16	17	18	19	20	21	22	23	24	25	26
	29	30	31	32	33	34	35	36	37	38	39	40
	43	44	45	46	47	48	49	50	51	52		

<b>Annually</b>	<b>Drain process fluid and replace with fresh fluid.</b>	
	<b>Check for fluid leaks throughout product and application.</b>	
	<b>Check filters for accumulation of particulate matter.</b>	
	<b>Clear any dust and debris from inside the product.</b>	
	<b>Vacuum out electrical box</b>	



# Annex I-1

## DOCUMENT DETAILS

Date	30/AUG/2021	Author(s)	MJH	Page	1 / 1	Revision	3
------	-------------	-----------	-----	------	-------	----------	---

## GENERIC MAINTENANCE FOR TECHNICIANS



Warning; Opening the refrigeration system may expose the operative to toxic and corrosive compounds (HF). Take protective measures including suitable eye protection.



Warning; Gases may exceed 300 psi (20 bar) during operation.



Warning; Refrigerants do not support combustion (A1-class), but do displace air (oxygen), presenting asphyxiation risk.



Warning; After switching off, the condenser cooling fan blades continue to rotate. Do not attempt servicing whilst the blades are rotating.



Warning; All chillers contain water and electricity in close proximity. Ensure the unit is isolated before service. This product is protected from overcurrent by fuses (or MCB) on the mains inlet. Never bypass the overcurrent protection.

Following service or repair by a trained technician;

- a) Ensure any electrical connections that may have been disturbed are given the 'tug-test'
- b) Ensure earth bonding conductors are re-attached.
- c) Ensure the correct fuses are in place.
- d) Ensure the mains cord being used is to specification, and is free from damage
- e) Subject the unit to a PAT test to ensure the unit is safe before running.
- f) Ensure there are no leaks inside or outside the unit.
- g) Using the wiring schematic for guidance, simulate faults to check each interlock's function.



Applied Thermal Control Ltd  
39 Hayhill Industrial Estate  
Barrow-upon-Soar, Loughborough  
LE12 8LD, United Kingdom  
+44 (0) 1530 839 998  
Service@thermalexchange.co.uk  
Support@app-therm.com

# Annex J-5

## DOCUMENT DETAILS

Date	03/FEB/2021	Author(s)	MJH	Page	1 / 1	Revision	01
------	-------------	-----------	-----	------	-------	----------	----

## CONFLICT MINERALS COMPLIANCE STATEMENT

Applied Thermal Control (ATC) adheres to and embraces the ethical values that support our everyday activities. As an expression of these principles and ethical values, ATC adheres to the principle of responsible sourcing of components containing precious and non-precious metals and metal salts in compliance with applicable laws and regulations.

The metals considered are Tantalum (Ta), Tungsten (W), Tin (Sn) and Gold (Au). ATC actively sources components from suppliers known to be reputable and could demonstrate compliance upon request with the Conflict Minerals acts and guidelines.

ATC uses Gold and Tin in electrical components, on PCBs and in rotating machinery, as governed by technical requirements of products. These metals could potentially originate from conflict mineral sites. As many of our suppliers do not purchase these metals direct from smelters, both they and ATC must rely heavily on information that will be provided by their suppliers to determine the source and chain of the metals in those products.

ATC is committed to working with its customers and supply chain to meet the customer's specification and requirements with regards to traceability, sourcing requirements and restrictions. ATC commits that, to the best of our knowledge, our suppliers are complying with the conflict minerals act as stated in their documentation. These statements are reviewed, and updates obtained as required.

Mitchell Howard, Technical Manager  
Signed in Coalville, UK, date 6/JUL/2020





Applied Thermal Control Ltd  
39 Hayhill Industrial Estate  
Barrow-upon-Soar, Loughborough  
LE12 8LD, United Kingdom  
+44 (0) 1530 839 998  
Service@thermalexchange.co.uk  
Support@app-therm.com

# Annex J-7

## DOCUMENT DETAILS

Date	03/FEB/2021	Author(s)	MJH	Page	1 / 1	Revision	01
------	-------------	-----------	-----	------	-------	----------	----

### WHAT IS THE REACH REGULATION 1907/2006?

REACH is a regulation of the European Union, adopted to improve the protection of human health and the environment from the risks that can be posed by chemicals, while enhancing the competitiveness of the EU chemicals industry. REACH places the burden of proof on companies. To comply with the regulation, companies must identify and manage the risks linked to the substances they manufacture and market in the EU. They have to demonstrate to ECHA how the substance can be safely used, and they must communicate the risk management measures to the users. If the risks cannot be managed, authorities can restrict the use of substances in different ways. In the long run, the most hazardous substances should be substituted with less dangerous ones. REACH stands for Registration, Evaluation, Authorization and Restriction of Chemicals. It entered into force on 1/JUN/2007.

### REACH 'ARTICLE' COMPLIANCE CONSIDERATIONS

#### REACH ANNEX XVII COMPLIANCE

Substances under Annex XVII are restricted either in full (not to be used at all) or for specific uses (can be used in some uses but cannot be used in identified uses).

Applied Thermal Control has contacted all our suppliers and to the best of our knowledge, none of the articles that we sell intentionally contain any of the Annex XVII substances currently on the candidate list in concentrations of >0.1% by weight.

#### REACH ANNEX XIV COMPLIANCE

Substances under Annex XIV require authorization to use in the EU after sunset date, require communication to downstream recipients when over threshold (0.1% w/w at article level) and require notification to ECHA when SVHC over threshold and imported over 1000kg annually and use not already registered.

Applied Thermal Control has contacted all our suppliers and to the best of our knowledge, none of the articles that we sell intentionally contain any of the Annex XVII substances currently on the candidate list in concentrations of >0.1% by weight.

#### SVHC LIST COMPLIANCE

Substances of Very High Concern (SVHC) require communication to downstream recipients when over threshold (0.1% w/w at the article level), notification to the European Chemicals Agency (ECHA) when SVHC over threshold and when imported over 1000kg annually and use not already registered.

Applied Thermal Control has contacted all our suppliers and to the best of our knowledge, none of the articles that we sell intentionally contain any of the Annex XVII substances currently on the candidate list in concentrations of >0.1% by weight.

### DECLARATION

Mitchell Howard, Technical Manager  
Signed in Barrow-upon-Soar, UK, date 15/JUL/2020



Applied Thermal Control Ltd  
39 Hayhill Industrial Estate  
Barrow-upon-Soar, Loughborough  
LE12 8LD, United Kingdom  
+44 (0) 1530 839 998  
Service@thermalexchange.co.uk  
Support@app-therm.com

# Annex J-8

## DOCUMENT DETAILS

Date	03/FEB/2021	Author(s)	MJH	Page	1 / 1	Revision	01
------	-------------	-----------	-----	------	-------	----------	----

### WHAT IS THE POPs REGULATION 2019/1021?

POPs stands for persistent organic pollutants. In Europe, the global Stockholm Convention is implemented through POPs legislation. POPs are organic substances that persist in the environment, accumulate in living organisms and pose a risk to our health and the environment. They can be transported by air, water or migratory species across international borders, reaching regions where they have never been produced or used. International risk management is necessary as no region can manage the risks posed by these substances alone.

The European Parliament (and Council) issued regulation 2019/1021 on 20/JUN/2019, and further amended (regulation 2020/784) on 8/APR/2020.

### POP<sub>s</sub> LISTED UNDER INITIAL REGULATION 2019/1021

**Pesticides;**

Aldrin, Chlordane, DDT, Dieldrin, Endrin, Heptachlor, Hexachlorobenzene, Mirex, Toxaphene.

**Industrial Chemicals;**

Hexachlorobenzene, Polychlorinated Biphenyls (PCBs).

**Industrial Chemical Byproducts;**

*Hexachlorobenzene byproducts;*

Polychlorinated dibenzo-p-dioxins and polychlorinated dibenzofurans (PCDD/PCDF), and PCBs.

### POP<sub>s</sub> LISTED UNDER AMENDMENT 2020/784

Perfluorooctanoic acid (PFOA), its salts and PFOA-related compounds.

### POP<sub>s</sub> COMPLIANCE STATEMENT

We certify that to the best of our knowledge, based upon up-to-date information from our suppliers, all products supplied by Applied Thermal Control are fully POPs compliant in accordance with regulations and amendments above mentioned.

### DECLARATION

Mitchell Howard, Technical Manager  
Signed in Barrow-upon-Soar, UK, date 27/AUG/2020



Applied Thermal Control Ltd  
39 Hayhill Industrial Estate  
Barrow-upon-Soar, Loughborough  
LE12 8LD, United Kingdom  
+44 (0) 1530 839 998  
Service@thermalexchange.co.uk  
Support@app-therm.com

## Operating Manual; Declarations & Approvals

# Annex J-10

### DOCUMENT DETAILS

Date	03/FEB/2021	Author(s)	MJH	Page	1 / 1	Revision	02
------	-------------	-----------	-----	------	-------	----------	----

### WHAT IS THE RoHS DIRECTIVE?

The RoHS Directive places restrictions on the use of certain hazardous substances in electrical and electronic equipment (EEE). RoHS compliance has been required for many years, however in 2014 it became a mandatory requirement under CE Marking. ATC products do not clearly fall within any of the existing categories of equipment, but as of 23/JUL/2019, all EEE not covered falls within scope of the directive. In contrast to RoHS 1, RoHS 2 is a CE marking Directive, and requires, for finished EEE, the use of the CE mark on the product to show compliance. The responsibility for affixing the CE mark resides with the manufacturer.

### RoHS 1 2002/95/EC

Adopted in February 2003 by the EU and taking effect on 1/JUL/2006, RoHS 1 restricted the use of 6 hazardous materials;

- 1) Lead (Pb)
- 2) Mercury (Hg)
- 3) Cadmium (Cd)
- 4) Hexavalent Chromium (Cr6+)
- 5) Polybrominated Biphenyls (PBB)
- 6) Polybrominated Diphenyl Ether (PBDE)

*We certify that to the best of our knowledge, based upon up-to-date information from our suppliers, all products supplied by Applied Thermal Control are fully RoHS 1 compliant.*

### RoHS 2 2011/65/EU

Adopted in July 2011 by the EU and taking effect on 2/JAN/2013, RoHS 2 expands the scope of RoHS 1 by adding new categories. RoHS 2 compliance is required to CE mark the product. Compliance with RoHS 2 is mandatory from 22/JUL/2019.

*We certify that to the best of our knowledge, based upon up-to-date information from our suppliers, all products supplied by Applied Thermal Control are fully RoHS 2 compliant.*

### RoHS 3 2015/863/EU

Adopted in 2015 by the EU and taking effect from 22/JUL/2019, RoHS 3 adds four additional substances to RoHS 1's list.

- 1) Bis(2-Ethylhexyl) phthalate (DEHP): < 1000 ppm
- 2) Benzyl butyl phthalate (BBP): < 1000 ppm
- 3) Dibutyl phthalate (DBP): < 1000 ppm
- 4) Di-isobutyl phthalate (DIBP): < 1000 ppm

*We certify that to the best of our knowledge, based upon up-to-date information from our suppliers, all products supplied by Applied Thermal Control are fully RoHS 3 compliant.*

### DECLARATION

Mitchell Howard, Technical Manager  
Signed in Barrow-upon-Soar, UK, date 11/NOV/2020



Applied Thermal Control Ltd  
 39 Hayhill Industrial Estate  
 Barrow-upon-Soar, Loughborough  
 LE12 8LD, United Kingdom  
 +44 (0) 1530 839 998  
 Service@thermalexchange.co.uk  
 Support@app-therm.com

Operating Manual; Declarations & Approvals

# Annex J-12

**DOCUMENT DETAILS**

Date	30/AUG/2021	Author(s)	MJH	Page	1 / 1	Revision	1
------	-------------	-----------	-----	------	-------	----------	---



**UKCA DECLARATION OF CONFORMITY (DoC)**

Demand created by; The Product Safety and Metrology etc. (Amendment etc.) (EU Exit) Regulations 2019

**REGISTERED BUSINESS ADDRESS**

Applied Thermal Control Ltd, 39 Hayhill Industrial Estate, Barrow-upon-Soar, Loughborough, LE12 8LD, UK.

**AUTHORISATION TO COMPILE THE TECHNICAL FILE**

Mitchell Howard, Applied Thermal Control Ltd, 39 Hayhill Industrial Estate, Barrow-upon-Soar, Loughborough, LE12 8LD, UK.

**DESCRIPTION & IDENTIFICATION OF MACHINERY**

Generic denomination;	R-Series
Function;	Recirculating chiller
Model;	All with 'R' prefix.
Type;	Air-cooled or water-cooled vapour compression-based.
Serial number;	
Commercial name;	As above.

**NOTIFIED BODY**

Not applicable

**QUALITY ASSURANCE SYSTEM**

QMS International Ltd, Muspole Court, Muspole Street, Norwich, NR3 1DJ, United Kingdom.  
 ASCB Registered; 201409-2

**DECLARATION**

The manufacturer declares that the machinery described above is in conformity with the relevant statutory requirements applicable to the specific product. The manufacturer takes full responsibility for the product's compliance.

- Supply of Machinery (Safety) Regulations 2008
- Electromagnetic Compatibility Regulations 2016
- Electrical Equipment (Safety) Regulations 2016
- The Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment Regulations 2012

**PERSON EMPOWERED TO DRAW UP DECLARATION**

Robert Poniatowski, CEO  
 Signed in Barrow-upon-Soar, UK, date 30/AUG/2021



Applied Thermal Control Ltd  
 39 Hayhill Industrial Estate  
 Barrow-upon-Soar, Loughborough  
 LE12 8LD, United Kingdom  
 +44 (0) 1530 839 998  
 Service@thermalexchange.co.uk  
 Support@app-therm.com

Operating Manual; Declarations & Approvals

# Annex J-13

**DOCUMENT DETAILS**

Date	30/AUG/2021	Author(s)	MJH	Page	1 / 1	Revision	1
------	-------------	-----------	-----	------	-------	----------	---

**EU DECLARATION OF CONFORMITY**

Document layout; Governed by Machinery Directive 2006/42/EC, Annex II.

**REGISTERED BUSINESS ADDRESS**

Applied Thermal Control Ltd, 39 Hayhill Industrial Estate, Barrow-upon-Soar, Loughborough, LE12 8LD, UK.

**AUTHORISATION TO COMPILE THE TECHNICAL FILE**

Mitchell Howard, Applied Thermal Control Ltd, 39 Hayhill Industrial Estate, Barrow-upon-Soar, Loughborough, LE12 8LD, UK.

**DESCRIPTION & IDENTIFICATION OF MACHINERY**

Generic denomination;	R-Series
Function;	Recirculating chiller
Model;	All with 'R' prefix.
Type;	Air-cooled or water-cooled vapour compression-based.
Serial number;	
Commercial name;	As above.

**NOTIFIED BODY**

Not applicable

**QUALITY ASSURANCE SYSTEM**

QMS International Ltd, Muspole Court, Muspole Street, Norwich, NR3 1DJ, United Kingdom.  
 ASCB Registered; 201409-2

**DECLARATION**

The manufacturer declares that the machinery described above fulfils all the relevant provisions of the;

- Machinery Directive 2006/42/EC.
- EMC Directive 2014/30/EU, via harmonised standards;
  - IEC 61000-6-2:2005 (Immunity for industrial environments).
  - IEC 61000-6-4:2006 +A1:2011 (Emission for industrial environments).
- Low Voltage Directive 2014/35/EU.
- RoHS Directive 2011/65/EU (*RoHS 2*);
  - The machinery above contains no Lead (Pb), Mercury (Hg), Cadmium (Cd), Hexavalent Chromium (Cr6+), Polybrominated Biphenyls (PBB) or Polybrominated Diphenyl Ether (PBDE).
- RoHS Directive (EU) 2015/863 (*RoHS 3*);
  - Bis(2-Ethylhexyl) phthalate (DEHP): < 1000 ppm
  - Benzyl butyl phthalate (BBP): < 1000 ppm
  - Dibutyl phthalate (DBP): < 1000 ppm
  - Diisobutyl phthalate (DIBP): < 1000 ppm

**PERSON EMPOWERED TO DRAW UP DECLARATION**

Robert Poniatowski, CEO

Signed in Barrow-upon-Soar, UK, date 30/AUG/2021



## DOCUMENT DETAILS

Date	03/FEB/2021	Author(s)	RW, MJH	Page	1 / 1	Revision	02
------	-------------	-----------	---------	------	-------	----------	----

## WARRANTY TERMS

Please visit the website warranty registration page to ensure ATC can offer you the best possible support;

**<https://www.app-therm.com/warranty-registration/>**

**a) For how long is my ATC product warrantied?**

ATC provides a comprehensive return to base 2-year parts, 1-year labor warranty from delivery as standard on all new equipment, provided it has been installed and operated in accordance with the manual.

**b) Where will ATC fulfill the product warranty?**

ATC's standard warranty terms are Return to Base (RTB) – issues with chillers are often easily solvable over the phone or email, or by reviewing ATC's technical guidance on the web and in the product manual. On occasion, at the discretion of ATC, goods may be serviced on site FOC or a service loan unit may be supplied. Warranty cover excludes the cost of travel by engineers and loan unit rental charges. Obtaining onsite service for a product, even in full warranty, is a chargeable service.

**c) Who is liable for shipping charges in the event of warranty failure?**

During the **first year** of the warranty period, freight costs for shipping to ATC are for the customer's account. Freight costs for shipping from ATC are for ATC's account.

During the **second year** of the warranty, freight costs to and from ATC are for the customer's account.

**d) I'm experiencing problems with my chiller. It's within warranty – what do I do next?**

Contact ATC to discuss the issue you are having. The contact details in the header of this document are an ideal place to start. Be sure to have your model number and serial number on-hand to aid those attempting to solve remotely.

**e) Telephone support couldn't fix my chiller – what do I do next?**

An RMA form must be completed. This allows both the end-user and ATC to clarify your details, to set the party responsible for shipping costs, and to set a different return address if desired. Shipping advice is provided, and the end-user must sign a declaration that states the unit is safe to handle. Return the form by email for fastest response.

**f) What happens if my chiller failed outside warranty or requires non-warranty repair work?**

A purchase order will be requested to cover an initial inspection – this will only be invoiced if the inspection shows there is no fault. If packaging is required, i.e. a crate, a separate charge will be levied. If the end user prefers ATC to arrange a collection, a shipping charge may be levied.

**g) Our process must continue running – can we have a loan unit whilst our chiller is in repair?**

ATC hold several standard air-cooled chillers at the factory for the sole purpose of offering for loan – these are available on a first-come, first-serve basis. Models up-to 3kW capacity are available.



# Annex M-9

## DOCUMENT DETAILS

Date	9/MAR/2022	Author(s)	MJH	Page	1 / 1	Revision	2
------	------------	-----------	-----	------	-------	----------	---

## RECOMMENDED SPARES FOR R10 & R20

Recommended spares include all rotating machinery (i.e. motors, fans), all sacrificial elements (i.e. fuses) and parts that users interact with (dials, fittings). Pricing is available from [service@thermalexchange.com](mailto:service@thermalexchange.com).

### R10

PN	Description	QTY
EA512	FUSE – T6.3A H250V UL-OK (FS1+2)	2
EA514	FUSE – T1.6A H250V UL-OK (FS3)	1
EA517	MOTOR – Standard 1/2~ P2=370W pump motor for P5 P10 P17	1
EA781	SWITCH – level switch threaded	1
63-500	FAN – Ø150mm, 1400RPM	2
WA679	COUPLING – bronze double flat male-to-male pump-to-motor	1
WA333	PUMP – P3 rotary vane	1
WA335	PUMP – P5 rotary vane	1
WA319	PUMP – P17 rotary vane	1
WA318	PUMP – P10 rotary vane	1
EA894	CONTROLLER – Ascon Tecnologic KR1	1
EA586	SENSOR – PT100	1
RA210	COMPRESSOR – KT1 / R10 1kW nominal duty (-0spec model)	1
EA898	RELAY – Solid state relay (SSR) DIN rail mount	1
EA496	RELAY – Power relay 30A	1
EA737	COIL – Solenoid valve coil	1
EA044	CORDSET – UK BS1363 to C13, 2m, right angle head (to suit appliance C14)	1

### R20

PN	Description	QTY
EA847	FUSE – T10A H250V UL-OK (FS1+2)	2
EA769	FUSE – T4A H250V UL-OK (FS3)	1
EA512	FUSE – T6.3A H250V UL-OK (FS4)	1
EA506	FUSE – T2A H250V UL-OK (FS5)	1
EA517	MOTOR – Standard 1/2~ P2=370W pump motor for P5 P10 P17	1
EA781	SWITCH – level switch threaded	1
63-500	FAN – Ø150mm, 1400RPM	4
WA679	COUPLING – bronze double flat male-to-male pump-to-motor	1
WA333	PUMP – P3 rotary vane	1
WA335	PUMP – P5 rotary vane	1
WA319	PUMP – P17 rotary vane	1
WA318	PUMP – P10 rotary vane	1
EA894	CONTROLLER – Ascon Tecnologic KR1	1
EA586	SENSOR – PT100	1
RA226	COMPRESSOR – K1 / R20 2kW nominal duty (-0spec model)	1
EA898	RELAY – Solid state relay (SSR) DIN rail mount	1
EA496	RELAY – Power relay 30A	1
EA737	COIL – Solenoid valve coil	1
EA044	CORDSET – UK BS1363 to C13, 2m, right angle head (to suit appliance C14)	1

## SAFETY DATA SHEET

according to Regulation (EU) 2015/830

Page 1/7

### Harp<sup>®</sup> 134a

Revision 0  
Revision date 2017-06-21

#### SECTION 1: Identification of the substance/mixture and of the company/undertaking

##### 1.1. Product identifier

<b>Product name</b>	Harp <sup>®</sup> 134a
<b>REACH Registration Number</b>	01-2119459374-33
<b>CAS No.</b>	811-97-2
<b>EC No.</b>	212-377-0

##### 1.2. Relevant identified uses of the substance or mixture and uses advised against

<b>Product Use</b>	[SU3] Industrial uses: Uses of substances as such or in preparations at industrial sites; -----
	[SU22] Professional uses: Public domain (administration, education, entertainment, services, craftsmen);
<b>Restricted use</b>	[SU21] Consumer uses: Private households (= general public = consumers);
<b>Description</b>	Gas.

##### 1.3. Details of the supplier of the safety data sheet

<b>Company</b>	Harp International Limited
<b>Address</b>	Gellihirion Industrial Estate Pontypridd Rhondda Cynon Taff CF37 5SX UK
<b>Web</b>	www.harpintl.com
<b>Telephone</b>	+44 (0)1443 842 255
<b>Fax</b>	+44 (0)1443 841 805
<b>Email</b>	harp@harpintl.com
<b>Email address of the competent person</b>	safety@harpintl.com

##### 1.4. Emergency telephone number


<b>Emergency telephone number</b>	+44 (0) 1270 502891 24 Hours
-----------------------------------	---------------------------------

#### SECTION 2: Hazards identification

##### 2.1. Classification of the substance or mixture

<b>2.1.2. Classification - EC 1272/2008</b>	Compressed gas: H280;
---	-----------------------

##### 2.2. Label elements

<b>Hazard pictograms</b>	
--------------------------	---



## Harp® 134a

Revision 0

Revision date 2017-06-21

## 2.2. Label elements

Signal Word	Warning
Hazard Statement	Compressed gas: H280 - Contains gas under pressure; may explode if heated.
Precautionary Statement: Storage	P410+P403 - Protect from sunlight. Store in a well-ventilated place.

## 2.3. Other hazards

Other hazards	Asphyxiant in high concentrations. May cause cold burns/frostbite.
---------------	--

## SECTION 3: Composition/information on ingredients

## 3.1. Substances

## 67/548/EEC / 1999/45/EC

Chemical Name	Index No.	CAS No.	EC No.	REACH Registration Number	Conc. (%w/w)	Classification	M-factor.
Harp® 134a (1,1,1,2-Tetrafluoroethane (HFC 134a))		811-97-2	212-377-0	01-2119459374-33	90 - 100%		

## EC 1272/2008

Chemical Name	Index No.	CAS No.	EC No.	REACH Registration Number	Conc. (%w/w)	Classification	M-factor.
Harp® 134a (1,1,1,2-Tetrafluoroethane (HFC 134a))		811-97-2	212-377-0	01-2119459374-33	90 - 100%	Compressed gas: H280;	

## SECTION 4: First aid measures

## 4.1. Description of first aid measures

Inhalation	Move the exposed person to fresh air.
Eye contact	Rinse immediately with plenty of water.
Skin contact	Frostbite: treat as thermal burns.
Ingestion	Ingestion is not considered a potential route of exposure.

## 4.2. Most important symptoms and effects, both acute and delayed

Inhalation	Seek medical attention if irritation or symptoms persist.
Eye contact	Seek medical attention if irritation or symptoms persist.
Skin contact	Frostbite: treat as thermal burns.
Ingestion	Ingestion is not considered a potential route of exposure.

## 4.3. Indication of any immediate medical attention and special treatment needed

Inhalation	If you feel unwell, seek medical advice (show the label where possible).
Eye contact	Seek medical attention if irritation or symptoms persist.
Skin contact	Seek medical attention if irritation or symptoms persist.
Ingestion	Ingestion is not considered a potential route of exposure.

## SECTION 5: Firefighting measures

## 5.1. Extinguishing media

	Use extinguishing media appropriate to the surrounding fire conditions.
--	---

## 5.2. Special hazards arising from the substance or mixture

	This product is not flammable in air under ambient conditions of temperature and pressure. Certain mixtures of the product and air under pressure may be flammable. At high temperature :, Thermal decomposition giving toxic and corrosive products :, Gaseous hydrogen fluoride (HF)., Carbon oxides.
--	---

## 5.3. Advice for firefighters

	Wear self contained breathing apparatus and protective clothing. Cool containers / tanks with water spray. Ensure a system for the rapid emptying of containers. In case of fire nearby, remove
--	---

## Harp® 134a

Revision 0

Revision date 2017-06-21

## 5.3. Advice for firefighters

exposed containers.

**SECTION 6: Accidental release measures**

## 6.1. Personal precautions, protective equipment and emergency procedures

Ensure adequate ventilation of the working area. Avoid contact with skin and eyes. Evacuate personnel to a safe area. Wear self contained breathing apparatus and protective clothing. Vapours are heavier than air.

## 6.2. Environmental precautions

Do not release into the environment.

## 6.3. Methods and material for containment and cleaning up

Recovery: Allow to evaporate.  
Elimination: See chapter 13.

## 6.4. Reference to other sections

See section 8. EXPOSURE CONTROLS / PERSONAL PROTECTION for further information.  
See section 13. DISPOSAL CONSIDERATIONS for further information.

**SECTION 7: Handling and storage**

## 7.1. Precautions for safe handling

Ensure adequate ventilation of the working area. Avoid contact with eyes and skin. Adopt best Manual Handling considerations when handling, carrying and dispensing. Keep away from sources of ignition - No smoking. Do not eat, drink or smoke in areas where this product is used or stored. When using do not eat or drink. Wash hands after handling the product.

## 7.2. Conditions for safe storage, including any incompatibilities

Keep containers tightly closed. Keep in a cool, dry, well ventilated area. Store in correctly labelled containers. Keep away from sources of ignition - No smoking. Store out of direct sunlight. Storage temperature: <45°C.

## Suitable packaging

Stainless steel. Steel.

## 7.3. Specific end use(s)

See section 1.2. Relevant identified uses of the substance or mixture and uses advised against for further information.

**SECTION 8: Exposure controls/personal protection**

## 8.1. Control parameters

Occupational exposure controls.

## 8.1.1. Exposure Limit Values

Harp® 134a (1,1,1,2-Tetrafluoroethane (HFC 134a))	WEL 8-hr limit ppm: 1000	WEL 8-hr limit mg/m3: 4240
	WEL 15 min limit ppm: -	WEL 15 min limit mg/m3: -
	WEL 8-hr limit mg/m3 total - inhalable dust:	WEL 15 min limit mg/m3 total - inhalable dust:
	WEL 8-hr limit mg/m3 total - respirable dust:	WEL 15 min limit mg/m3 total - respirable dust:


## 8.2. Exposure controls

## Harp® 134a

Revision 0

Revision date 2017-06-21

## 8.2. Exposure controls

	
8.2.1. Appropriate engineering controls	Ensure adequate ventilation of the working area.
8.2.2. Individual protection measures	Wear protective clothing.
Eye / face protection	Approved safety goggles.
Skin protection - Handprotection	Wear suitable gloves.
Skin protection - Other	Wear suitable protective clothing.
Respiratory protection	Wear suitable respiratory equipment when necessary.
Occupational exposure controls	Keep away from food, drink and animal feedingstuffs.

## SECTION 9: Physical and chemical properties

## 9.1. Information on basic physical and chemical properties

Appearance	Gas
Colour	Colourless
Odour	Slight
Odour threshold	No data available
Freezing Point	No data available
Evaporation rate	No data available
Water solubility	No data available
Fat Solubility	No data available
Soluble in	No data available
Partition coefficient (n-octanol/water)	No data available
Partition coefficient	No data available
Autoignition temperature	> 743 °C
Decomposition temperature	> 370 °C
Vapour pressure	= 0.574 MPa
Vapour density	= 4.24 kg/m <sup>3</sup>
Relative density	= 1.21 (H <sub>2</sub> O = 1 @ 20 °C)
Initial boiling point	- 26 °C
Melting point	- 108 °C
Flash point	Not applicable.
pH	Not applicable.
Flammability (solid, gas)	Not applicable.
Viscosity	Not applicable.
Explosive properties	Not applicable.
Oxidising properties	Not applicable.
Solubility	No data available

## 9.2. Other information

## Harp® 134a

Revision 0

Revision date 2017-06-21

## 9.2. Other information

VOC (Volatile organic compounds)	Not relevant
Conductivity	No data available
Surface tension	No data available
Gas group	No data available
Benzene Content	No data available
Lead content	No data available

## SECTION 10: Stability and reactivity

## 10.1. Reactivity

	Stable under normal conditions. The gaseous product in presence of air can form, under certain conditions of temperature and pressure, a flammable mixture.
--	---

## 10.2. Chemical stability

	Stable under normal conditions. The gaseous product in presence of air can form, under certain conditions of temperature and pressure, a flammable mixture.
--	---

## 10.3. Possibility of hazardous reactions

	No data is available on this product.
--	---------------------------------------

## 10.4. Conditions to avoid

	Keep away from heat and sources of ignition. Avoid contact with flames and red hot metallic surfaces.
--	---

## 10.5. Incompatible materials

	Alkaline hydroxides. Alkaline earth metals. Strong oxidising agents. Finely divided metals.
--	---

## 10.6. Hazardous decomposition products

	At high temperature ; Thermal decomposition giving toxic and corrosive products ; Gaseous hydrogen fluoride (HF)., Carbon oxides. Decomposition temperature: >370°C.
--	--

## SECTION 11: Toxicological information

## 11.1. Information on toxicological effects

Acute toxicity	Slightly harmful by inhalation. As with other volatile aliphatic halogenated compounds, through vapour accumulation and/or inhalation of large quantities, the product can cause ; Loss of consciousness and cardiac disorders aggravated by stress and lack of oxygen, risk of mortality.
Skin corrosion/irritation	Ejection of liquefied gas : frostbite possible.
Serious eye damage/irritation	Ejection of liquefied gas : frostbite possible.
Respiratory or skin sensitisation	No data available.
Germ cell mutagenicity	No data available.
Carcinogenicity	No data available.
Reproductive toxicity	No data available.
STOT-single exposure	No data available.
STOT-repeated exposure	No data available.
Aspiration hazard	No data available.
Repeated or prolonged exposure	No data available.

## 11.1.4. Toxicological Information

	No data available
--	-------------------

## SECTION 12: Ecological information

## Harp® 134a

Revision 0  
Revision date 2017-06-21

## 12.1. Toxicity

No data available

## 12.2. Persistence and degradability

Not readily biodegradable.

## 12.3. Bioaccumulative potential

Does not bioaccumulate.

## Partition coefficient

Harp® 134a No data available

## 12.4. Mobility in soil

No data is available on this product.

## 12.5. Results of PBT and vPvB assessment

This substance is not considered to be persistent, bioaccumulating, toxic (PBT), nor very persistent, very bioaccumulating (vPvB).

## 12.6. Other adverse effects

Global warming potential (GWP): Global warming potential with respect to CO<sub>2</sub> = 1430 (IPCC Assessment Report 4).  
Ozone depletion potential: Ozone depletion potential; ODP; (R-11 = 1), Value:.

## SECTION 13: Disposal considerations

## 13.1. Waste treatment methods

Dispose of in compliance with all local and national regulations.

## Disposal methods

Contact a licensed waste disposal company.

## SECTION 14: Transport information

## Hazard pictograms



## 14.1. UN number

UN3159

## 14.2. UN proper shipping name

1,1,1,2-TETRAFLUOROETHANE

## 14.3. Transport hazard class(es)

ADR/RID	2
Subsidiary risk	-
IMDG	2.2
Subsidiary risk	-
IATA	2.2
Subsidiary risk	-

## 14.4. Packing group

Packing group -

## Harp® 134a

Revision 0

Revision date 2017-06-21

## 14.5. Environmental hazards

Environmental hazards	No
Marine pollutant	No

## 14.6. Special precautions for user

	No data is available on this product.
--	---------------------------------------

## 14.7. Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code

	No data is available on this product.
--	---------------------------------------

## ADR/RID

Hazard ID	20
Tunnel Category	(C/E)

## IMDG

EmS Code	F-C S-V
----------	---------

## IATA

Packing Instruction (Cargo)	200
Maximum quantity	150 kg
Packing Instruction (Passenger)	200
Maximum quantity	75 kg

## SECTION 15: Regulatory information

## 15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

Regulations	REGULATION (EC) No 1907/2006 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 18 December 2006 concerning the Registration, Evaluation, Authorisation and Restriction of Chemicals (REACH), establishing a European Chemicals Agency, amending Directive 1999/45/EC and repealing Council Regulation (EEC) No 793/93 and Commission Regulation (EC) No 1488/94 as well as Council Directive 76/769/EEC and Commission Directives 91/155/EEC, 93/67/EEC, 93/105/EC and 2000/21/EC.
-------------	---

## 15.2. Chemical safety assessment

	No data is available on this product.
--	---------------------------------------

## SECTION 16: Other information

## Other information

Text of Hazard Statements in Section 3	Compressed gas: H280 - Contains gas under pressure; may explode if heated.
--	--

## Further information

	The information supplied in this Safety Data Sheet is designed only as guidance for the safe use, storage and handling of the product. This information is correct to the best of our knowledge and belief at the date of publication however no guarantee is made to its accuracy. This information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any other process.
--	--

**HARP® R407C**

Version: CLP01

Date: July 2012

Page 1 of 8

**1. Identification of the substance / preparation and company / undertaking**

Product name	R407C	
REACH registration numbers	1,1,1,2-Tetrafluoroethane Pentafluoroethane Difluoromethane	01-2119459374-33 01-2119485636-25 Deadline not yet expired
Company	Harp International Ltd Gellihirion Industrial Estate Pontypridd Rhondda Cynon Taff CF37 5SX Tel: +44 (0) 1443 842255 Fax: +44 (0) 1443 841805 Email: harp@harpintl.com	
Emergency phone number	+44 (0) 1270 502891 (24 hour)	
Use	Refrigeration	

**2. Hazards identification**

**EC Classification**

EC Directive 67/548/EEC or 1999/45/EC	Not classified as hazardous
Regulation (EC) No. 1272/2008 (CLP)	Gases under pressure – Liquefied gas

**Label Elements**

Name on label	
Hazardous components	1,1,1,2-Tetrafluoroethane (R134a) Pentafluoroethane (R125) Difluoromethane (R32)
Hazard statement(s)	H280: Contains gas under pressure; may explode if heated
Signal word(s)	Warning
Hazard pictogram(s)	



Precautionary statement(s)	
Storage	P410 + P403: Protect from sunlight. Store in a well-ventilated place.

# SAFETY DATA SHEET

According to Regulation (EC) No.1907/2006

  
HARP  
INTERNATIONAL

## HARP® R407C

Version: CLP01

Date: July 2012

Page 2 of 8

### 3. Composition / information on ingredients

#### Concentration

Substance name	CAS No.	EC No.	Concentration
1,1,1,2-Tetrafluoroethane (R134a)	811-97-2	212-377-0	ca. 52%
Pentafluoroethane (R125)	354-33-6	206-557-8	ca. 25%
Difluoromethane (R32)	75-10-5	200-839-4	ca. 23%

#### Hazardous components according to Regulation (EC) 1272/2008 as amended

Substance name	Hazard class	Hazard category	H Phrases
1,1,1,2-Tetrafluoroethane (R134a)	Gases under pressure	Liquefied gas	H280
Pentafluoroethane (R125)	Gases under pressure	Liquefied gas	H280
Difluoromethane (R32)	Flammable gases	Category 1	H220
	Gases under pressure	Liquefied gas	H280

#### Hazardous components according to European Directive 67/548/EEC or 1999/45/EC as amended

Substance name	Classification	Hazard category	R-phrase(s)
Difluoromethane (R32)	F+	Extremely flammable	R12

### 4. First aid measures

<b>Inhalation</b>	Remove to fresh air. Oxygen or artificial respiration if needed. If symptoms persist, call a physician.
<b>Skin contact</b>	Allow to evaporate. Wash off with warm water. If symptoms persist, call a physician.
<b>Eye contact</b>	Immediately irrigate with eyewash solution or clean water, holding the eyelids apart for at least 10 minutes. Obtain immediate medical attention.
<b>Ingestion</b>	Unlikely route of exposure.
<b>Most important symptoms/effects, acute and delayed</b>	
<b>Inhalation</b>	In case of higher concentrations: narcosis, asphyxia, may cause cardiac arrhythmia.
<b>Skin contact</b>	Contact with liquid or refrigerated gas can cause cold burns and frostbite. Prolonged skin contact may defat the skin and produce dermatitis.
<b>Eye contact</b>	Causes frostbite burns to eyes. Symptoms: Lachrymation, redness, swelling of tissue, frostbite, burn.
<b>Ingestion</b>	Gas. Not applicable.

### 5. Fire-fighting measures

#### Extinguishing media

Suitable extinguishing media

As appropriate for surrounding fire. Keep fire exposed containers cool by spraying with water.

Unsuitable extinguishing media

None.



# SAFETY DATA SHEET

According to Regulation (EC) No.1907/2006

  
HARP  
INTERNATIONAL

## HARP® R407C

Version: CLP01

Date: July 2012

Page 3 of 8

### Specific hazards arising from the Chemical

The product is not flammable.  
Hazardous decomposition products formed under fire conditions.

### Special protective actions for Fire-Fighters

Wear self-contained breathing apparatus and protective suit  
Wear chemical resistant oversuit  
Special protective actions for fire-fighters  
In case of fire, use water spray  
Keep product and empty container away from heat and sources of ignition

## 6. Accidental release measures

### Personal precautions, protective equipment and emergency procedures

#### Advice for non-emergency personnel

Prevent further leakage or spillage if safe to do so  
Keep away from incompatible products

#### Advice for emergency responders

Immediately evacuate personnel to safe areas  
Keep people away from and upwind of spill/leak  
Wear self-contained breathing apparatus and protective suit  
Vapours are heavier than air and can cause suffocation by reducing oxygen available for breathing  
Suppress (knock down) gases/vapours/mists with a water spray jet  
Avoid spraying the leak source  
Ventilate area

#### Environmental precautions

Discharge into the environment must be avoided  
Inform the responsible authorities in case of gas leakage or of entry into waterways, soil or drains

#### Methods and materials for containment and cleaning up

Allow to evaporate  
Prevent product from entering drains

#### Reference to other sections

Refer to protective measures listed in sections 7 and 8.

## 7. Handling and storage

### Precautions for safe handling

Use only in well-ventilated areas  
Use only clean and dry utensils  
Keep away from water  
Preferably transfer by pump or gravity  
Keep away from incompatible products

### Conditions for storage, including incompatibilities

#### Storage

Keep only in the original container  
Store in a receptacle equipped with a vent  
Keep containers tightly closed in a cool, well-ventilated place  
Keep in properly labelled containers  
Keep in a bonded area  
Keep away from heat/sparks/open flames/hot surfaces. No smoking.  
Keep away from incompatible products

#### Packing material

Suitable material – steel cylinder

#### Specific use(s)

For further information, please contact supplier.

**HARP® R407C**

Version: CLP01

Date: July 2012

Page 4 of 8

**8. Exposure controls / personal protection**

**Control parameters**

Exposure limit values

Substance	Harp acceptable exposure limit	EH40 workplace exposure limits
1,1,1,2-Tetrafluoroethane	TWA = 1000 ppm	TWA = 1000 ppm / 4240 mg/m <sup>3</sup>
Pentafluoroethane	TWA = 1000 ppm	Not listed
Difluoromethane	TWA = 1000 ppm	Not listed

**Exposure controls**

Appropriate engineering controls	Ensure adequate ventilation Apply technical measures to comply with the occupational exposure limits
Respiratory protection	Self-contained breathing apparatus (EN 133) Wear self-contained breathing apparatus in confined spaces, in cases where the oxygen level is depleted, or in case of significant emissions Use only respiratory protection that conforms to international / national standards
Hand protection	Take note of the information given by the producer concerning permeability and break through times and of special workplace conditions (mechanical strain, duration of contact). Protective gloves Suitable material: Fluoroelastomer
Eye protection	Tightly fitted safety goggles
Skin and body protection	Wear suitable protective clothing If splashes are likely to occur, wear: apron, boots, Neoprene
Hygiene measures	Eye wash bottles or eye wash stations in compliance with applicable standards When using do not eat, drink or smoke Gloves, overalls and boots have to be double layered (protection against cold temperature). Handle in accordance with good industrial hygiene and safety practice
Environmental exposure controls	Dispose of rinse water in accordance with local and national regulations.

**9. Physical and chemical properties**

<b>Form</b>	Compressed liquefied gas
<b>Colour</b>	Colourless
<b>Odour</b>	Ether-like
<b>pH</b>	Neutral
<b>pKa</b>	Not applicable
<b>Melting point/freezing point</b>	-103°C (Pentafluoroethane)
<b>Boiling point/boiling range</b>	-44 to -37°C
<b>Flash point</b>	Not applicable
<b>Evaporation rate</b>	No data
<b>Flammability (solid, gas)</b>	The product is not flammable
<b>Flammability</b>	Not applicable

# SAFETY DATA SHEET

According to Regulation (EC) No.1907/2006

## HARP® R407C

Version: CLP01

Date: July 2012

Page 5 of 8

<b>Explosive properties</b>	Not explosive
<b>Vapour pressure</b>	10.35 bar at 20°C 21.94 bar at 50°C (Pentafluoroethane)
<b>Vapour density</b>	3.45
<b>Density</b>	Not applicable
<b>Relative density</b>	1.17 at 20°C
<b>Bulk density</b>	Not applicable
<b>Solubility</b>	430 mg/l at 25°C, water (Pentafluoroethane)
<b>Solubility/qualitative</b>	No data available
<b>Partition coefficient: n-octanol/water</b>	log Pow: 1.48, 20°C (pentafluoroethane)
<b>Auto-ignition temperature</b>	No data available
<b>Decomposition temperature</b>	No data
<b>Viscosity</b>	Not applicable
<b>Oxidizing properties</b>	Non oxidizer

### 10. Stability and reactivity

<b>Reactivity</b>	Risk of violent reaction
<b>Chemical stability</b>	Stable under recommended storage conditions
<b>Possibility of hazardous reactions</b>	Strong oxidizers, alkali metals and alkaline earth metals may cause fires or explosions.
<b>Conditions to avoid</b>	Heat
<b>Materials to avoid</b>	Light and/or alkaline metals, powdered metals, alkaline earth metals, oxidising agents
<b>Hazardous decomposition products</b>	Gaseous hydrogen fluoride (HF), Fluorophosgene The release of other hazardous decomposition products is possible

### 11. Toxicological information

<b>Acute toxicity</b>	
Acute oral toxicity	Not applicable
Acute inhalation toxicity	LC50, 4 h, >2,080,000 mg/m <sup>3</sup> (1,1,1,2-Tetrafluoroethane)
Acute dermal toxicity	Not applicable
<b>Skin corrosion/irritation</b>	Not applicable
<b>Serious eye damage/eye irritation</b>	Not applicable
<b>Respiratory or skin sensitization</b>	Guinea pig, did not cause sensitization on laboratory animals
<b>Mutagenicity</b>	In vitro tests did not show mutagenic effects (Pentafluoroethane) In vivo tests did not show mutagenic effects (Pentafluoroethane)
<b>Carcinogenicity</b>	Negative (1,1,1,2-Tetrafluoroethane)
<b>Toxicity for reproduction</b>	Developmental toxicity, rat, no observed effect (1,1,1,2-Tetrafluoroethane)
<b>Repeated dose toxicity</b>	Inhalation, after a single exposure, dog, cardiac sensitization following adrenergic stimulation Inhalation, rat, >=50000ppm, NOAEL (1,1,1,2-Tetrafluoroethane) Inhalation, repeated exposure, rat, >=50000ppm, NOAEL (Pentafluoroethane) Inhalation, 90-day, rat, 108 mg/m <sup>3</sup> , NOAEL (Difluoromethane)
<b>Other information</b>	No data available

## 12. Ecological information

### Toxicity

Fishes	Brachydanio rerio	LC50	96 h	>200 mg/l	1,1,1,3,3-pentafluorobutane
Fishes	Brachydanio rerio	LC50	96 h	Ca. 200 mg/l	1,1,1,3,3-pentafluorobutane
Crustaceans	Daphnia magna	EC50	48 h	>200 mg/l	1,1,1,3,3-pentafluorobutane
Crustaceans	Daphnia magna	NOEC	48 h	200 mg/l	1,1,1,3,3-pentafluorobutane
Algae	Selenastrum capricornutum	NOEC	72 h	13.2 mg/l	1,1,1,3,3-pentafluorobutane
Algae	Selenastrum capricornutum	EC50	72 h	>114 mg/l	1,1,1,3,3-pentafluorobutane
Terrestrial plants		NOEC	growth	>=6 g/m <sup>3</sup>	

### Persistence and degradability

#### Abiotic degradation

Air, indirect photo-oxidation. T ½ from 4.16 – 28.2 y

Conditions: sensitizer: OH radicals.

Degradation products: carbon dioxide (CO<sub>2</sub>) / hydrofluoric acid / TFA

#### Biodegradation

Aerobic, tested according to closed bottle test, chemical degradation, 2-5% after 28 d. Result: not readily biodegradable

Aerobic, tested according to biodegradation by methane oxidation. Result: not readily biodegradable (1,1,1,2-Tetrafluoroethane)

### Bioaccumulative potential

Bioaccumulative potential: log Pow 0.21-1.48. Result: does not bioaccumulate

### Mobility

Soil/sediments, adsorption, log KOC: from 1.05 – 1.7. Conditions: calculated value

Air, Henry's law constant (H), from 19.7 – 150 hPa.m<sup>3</sup>/mol, 20°C. Conditions: calculated value, considerable volatility

### Other adverse effects

Ozone depletion potential = 0

Result = no effect on stratospheric ozone

Ozone depletion potential; ODP; (R11 = 1)

Global Warming Potential = 0.25

Halocarbon global warming potential; HGWP; (R11 = 1) (1,1,1,2-Tetrafluoroethane)

## 13. Disposal considerations

### Waste disposal methods

In accordance with local and national regulations

Refer to manufacturer/supplier for information on recovery/recycling

### Contaminated packaging

To avoid treatments, as far as possible, use dedicated containers

Where possible recycling is preferred to disposal or incineration

## 14. Transport information

### International transport regulations

#### IATA-DGR

UN number

UN 3340

Class

2.2

ICAO-Labels

2.2 - Non-flammable, non-toxic gas

Proper shipping name

REFRIGERANT GAS R407C

# SAFETY DATA SHEET

According to Regulation (EC) No.1907/2006

## HARP® R407C

Version: CLP01

Date: July 2012

Page 7 of 8

### IMDG

UN number UN 3340  
Class 2.2  
IMDG-Labels 2.2 - Non-flammable, non-toxic gas  
HI/UN No. 3340  
EmS F-C, S-V  
Proper shipping name REFRIGERANT GAS R407C

### ADR

UN number UN 3340  
Class 2  
ADR/RID Labels 2.2 - Non-flammable, non-toxic gas  
HI/UN No. 20 / 3340  
Proper shipping name REFRIGERANT GAS R407C

### RID

UN number UN 3340  
Class 2  
ADR/RID Labels 2.2 - Non-flammable, non-toxic gas  
HI/UN No. 20 / 3340  
Proper shipping name REFRIGERANT GAS R407C

### ADN

UN number UN 3340  
Class 2  
ADR/RID Labels 2.2 – Non-flammable, non-toxic gas  
Proper shipping name REFRIGERANT GAS R407C

## 15. Regulatory information

### Applicable Laws or Regulations

- Regulation (EC) No 1907/2006 of the European Parliament and of the Council of 18 December concerning the Registration, Evaluation, Authorisation and Restriction of Chemicals (REACH) as amended
- Directive 1999/45/EC of the European Parliament and of the Council of 31 May 1999 concerning the approximation of laws, regulations and administrative provisions of the Member States relating to the classification, packaging and labelling of dangerous preparations, as amended
- Regulation (EC) No 1272/2008 of the European Parliament and of the Council of 16 December 2008 on classification, labelling and packaging of substances and mixtures, as amended
- Regulation (EC) No 166/2006 of the European Parliament and of the Council of 18 January 2006 concerning the establishment of a European Pollutant Release and Transfer Register and amending Council Directives 91/689/EEC and 96/61/EC
- Directive 2008/98/EC of the European Parliament and of the Council of 19 November 2008 on waste
- EH40/2005 Workplace Exposure Limits, as amended through 1, 10, 2007 (WEL's) published by the Health and Safety Executive (HSE). Issued under the Control of Substances Hazardous to Health Regulations, as amended

### Notification status

Inventory information	Status
Australian Inventory of Chemical Substances (AICS)	In compliance with inventory
Canadian Domestic Substances List (DSL)	In compliance with inventory
Inventory of Existing Chemical Substances (China) (IECS)	In compliance with inventory
Japanese Existing and New Chemical Substances (MITI List) (ENCS)	In compliance with inventory
New Zealand Inventory of Chemicals (NZIOC)	In compliance with inventory
Toxic Substance Control Act List (TSCA)	In compliance with inventory
EU List of Existing Chemical Substances (EINECS)	In compliance with inventory
Korean Existing Chemicals Inventory (KECI (KR))	In compliance with inventory
Philippine Inventory of Chemicals and Chemical Substances (PICCS)	In compliance with inventory

# SAFETY DATA SHEET

According to Regulation (EC) No.1907/2006

 HARP  
INTERNATIONAL

## HARP® R407C

Version: CLP01

Date: July 2012

Page 8 of 8

### 16. Other information

#### Full text of H-Statements referred to under section 3

H220	Extremely flammable gas
H280	Contains gas under pressure; may explode if heated

This datasheet was prepared in accordance with Regulation (EC) No. 1907/2006.

Information in this publication is believed to be accurate and is given in good faith but it is for the user to satisfy itself of the suitability for its own particular purpose. Accordingly, Harp International Limited gives no warranty as to the fitness of the product for any particular purpose and any implied warranty or condition, statutory or otherwise, is excluded except to the extent that such exclusion is prevented by law. Freedom under Patent, Copyright and Designs cannot be assumed. HARP® is a trademark, the property of Harp International Ltd.

Harp International Limited is registered in England & Wales No. 02432294. Registered office: Gellihirion Industrial Estate, Pontypridd, Rhondda Cynon Taff, CF37 5SX.

# SAFETY DATA SHEET

## HEXID A4 HEAT TRANSFER FLUID

Conforming to Directive 1907/2006/EC

### SECTION 1: IDENTIFICATION OF THE SUBSTANCE/MIXTURE AND OF THE COMPANY/UNDERTAKING

<b>1.1. Product Name</b>	Hexid A4
<b>1.2. Supplier</b>	Applied Thermal Control Limited 39 Hayhill Industrial Estate, Barrow upon Soar, Leicestershire, LE12 8LD. United Kingdom. www.app-therm.com
<b>1.3. Telephone Number</b>	+44(0)1530 839998
<b>1.4. Email</b>	<a href="mailto:sales@app-therm.com">sales@app-therm.com</a>
<b>1.5. Emergency Telephone Number</b>	+44(0)1530 839998
<b>1.6. Intended/Recommended Use</b>	Heat Transfer Fluid

### SECTION 2: HAZARDS IDENTIFICATION

- 2.1. Classification of the substance or mixture**  
The product is not classified as dangerous according to Regulation (EC) No. 1272/2008.  
This mixture is not classified as dangerous according to Directive 1999/45/EC.

### SECTION 3: COMPOSITION/INFORMATION ON INGREDIENTS

- 3.1. Chemical Nature** Water (CAS 7732-18-5), not classified.  
Propylene glycol (CAS 57-55-6) (REACH 01-2119456809-23)  
(EINECS 200-338-0) not classified.  
Fluorescein (trace) and biocide (trace) not classified.
- 3.2. Food Grade**

### SECTION 4: FIRST AID MEASURES

- General advise** No special precautions required. Treat symptomatically.
- 4.1. Eye Contact** Rinse thoroughly with plenty of water, also under the eyelids. Remove contact lenses after a few minutes and continue rinsing. If symptoms persist, call a physician.
- 4.2. Skin Contact** Wash off immediately with plenty of water. If skin irritation persists, call a physician.
- 4.3. Inhalation** Remove to fresh air. If symptoms persist, call a physician.
- 4.4. Ingestion** Rinse mouth with water. Never give anything by mouth to an unconscious person. If symptoms persist, call a physician.

### SECTION 5: FIREFIGHTING MEASURES

- 5.1. Extinguishing media**  
Use extinguishing measures that are appropriate to local circumstances and the surrounding environment. Water spray, foam, dry powder or CO<sub>2</sub>. Alcohol-resistant foam
- 5.2. Unsuitable extinguishing Media**  
High volume water jet. Do not use a solid water stream as it may scatter and spread fire.
- 5.3. Specific hazards during firefighting**  
In fire conditions, toxic decomposition products may be formed (see also section 10). In combustion, emits fumes, smoke, carbon dioxide (CO<sub>2</sub>) and carbon monoxide (CO). Heating will cause a pressure rise - with severe risk of bursting and explosion, Violent steam generation or eruption may occur upon application of direct water to hot liquids.
- 5.4. Advice for firefighters**  
In the event of fire, wear self-contained breathing apparatus. Wear personal protective equipment. Promptly isolate the scene by removing all persons from the vicinity of the incident if there is a fire. Keep containers cool by spraying with water if exposed to fire. Collect contaminated fire extinguishing water separately. This must not be discharged into drains. Burning fluids may be extinguished by dilution with water

# SAFETY DATA SHEET

## HEXID A4 HEAT TRANSFER FLUID

Conforming to Directive 1907/2006/EC

### SECTION 6: ACCIDENTAL RELEASE MEASURES

#### 6.1. Personal precautions

Use personal protective equipment. Avoid contact with skin and eyes. Keep unnecessary and unprotected personnel from entering the area.

#### 6.2. Precaution to protect the environment

Do not flush into surface water or sanitary sewer system. Avoid subsoil penetration.

#### 6.3. Clean-up procedures

Contain the spillage, soak up with non-combustible absorbent material, (e.g. sand, earth, diatomaceous earth, vermiculite) and transfer to a container for disposal according to local / national regulations (see section 13). Keep in suitable, closed containers for disposal. Dike the area of spill to prevent spreading and pump liquid to salvage tank. Treat recovered material as described in section 13 Disposal considerations.

### SECTION 7: HANDLING AND STORAGE

#### 7.1. Precautions for safe handling

Keep container tightly closed. Handle in accordance with good industrial hygiene and safety practice. Spills of these organic materials on hot fibrous insulations may lead to lowering of the auto-ignition temperatures possibly resulting in spontaneous combustion.

#### 7.2. Conditions for safe storage

Keep only in the original container.

### SECTION 8: EXPOSURE CONTROLS/PERSONAL PROTECTION

#### 8.1. Control parameters

Component: Propane-1,2-diol CAS-No. 57-55-6

Other Occupational Exposure Limit Values EH40 WEL, Time Weighted Average (TWA):, Total vapour and particulates.150 ppm, 474 mg/m<sup>3</sup>

EH40 WEL, Time Weighted Average (TWA):, Particulate.10 mg/m<sup>3</sup>

ELV (IE), Time Weighted Average (TWA):, Total vapour and particulates.150 ppm, 470 mg/m<sup>3</sup>

ELV (IE), Time Weighted Average (TWA):, Particulate.10 mg/m<sup>3</sup>

#### 8.2. Exposure controls/Appropriate engineering controls

Local exhaust. If this product contains ingredients with exposure limits, use process enclosures, local exhaust ventilation or other engineering controls to keep worker exposure below any recommended or statutory limits.

#### Personal protective equipment

Respiratory protection Suitable respiratory protective device Combination filter: A-P2

Filter Type Combined particulates and organic vapour type

Hand protection Category short time exposure Break through time > 10 min

Protective index Class 1 When prolonged exposure is expected: Break through time > 120 min

Protective index Class 4 Observe the information of the glove manufacturers on permeability.

Protective gloves should be chosen according to Workplace Safety Assessment.

Gloves recommended according to EN 374 (protection against chemicals).

Material Chemical resistant gloves made of butyl rubber or nitrile rubber category III according to EN 374.

### SECTION 9: PHYSICAL AND CHEMICAL PROPERTIES

9.1	Appearance at 20°C	Fluorescent green clear liquid
9.2	Odour	Almost odourless
9.3	Flash point	Boils without flashing
9.4	Ignition temperature	Not Available
9.5	Flammability Limit	Not Available
9.6	Oxidizing Properties	Not Available
9.7	Auto flammability	450°C
9.8	Density at 25°C	~1.036g/cm <sup>3</sup>
9.9	pH (as is)	7
9.10	Boiling point	102°C
9.7	Auto flammability	450°C
9.8	Solubility in water	Miscible
9.9	Freezing point	-21°C



# SAFETY DATA SHEET

## HEXID A4 HEAT TRANSFER FLUID

Conforming to Directive 1907/2006/EC

9.10	Specific Heat Capacity	3.78kJ/kg °K
9.11	Viscosity, Kinetic, at 25°C	3.51mPa.s

---

### SECTION 10: STABILITY AND REACTIVITY

---

#### 10.1. Reactivity

Stable under recommended storage conditions. No dangerous reaction known under conditions of normal use.

#### 10.2. Chemical stability

No decomposition if stored and applied as directed. Stable under recommended storage conditions. Hygroscopic.

#### 10.3. Hazardous reactions

Hazardous polymerisation does not occur.

#### 10.4. Conditions to avoid

Generation of gas from decomposition causes pressure in closed systems. Keep away from direct sunlight. Avoid high temperatures. Avoid temperatures exceeding the decomposition temperature. Avoid UV light.

#### 10.5. Materials to avoid

Strong acids, Strong bases, Strong oxidizing agents.

#### 10.6. Hazardous decomposition products

Aldehydes, Alcohols, Ether, Organic acids.

---

### SECTION 11: TOXICOLOGICAL INFORMATION

---

#### 11.1. Toxicity Oral

LD50 : > 20000 mg/kg (rat) This product can present a small hazard if large quantities are swallowed.

#### 11.2. Inhalation

LC50 : 6.15 mg/l (rat; 4 h; vapour) At ambient temperature the exposure to vapours is minimal due to a low volatility rate. Inhalation may cause irritation to the nose, throat, upper respiratory tract and lungs. No deaths occurred

#### 11.3. Dermal

LD50 : > 20000 mg/kg (rabbit) Prolonged skin contact is unlikely to result in absorption of harmful amounts. Skin irritation by prolonged exposure is unlikely. Repeated contact may cause flaking and softening of skin.

#### 11.4. Eyes

Slight irritation is possible. Direct contact with eyes may cause temporary irritation. Corneal injury is unlikely.

#### 11.5. Sensitisation

Patch test on human volunteers did not demonstrate sensitisation properties.

#### 11.6. CMR Carcinogenicity

Animal testing did not show any carcinogenic effects. Information given is based on data obtained from similar substances.

#### 11.7. Mutagenicity

No data available.

#### 11.8. Reproductive toxicity

No data available.

#### 11.9. Specific Target Organ Toxicity

Single exposure no data available. Repeated exposure no data available.

#### 11.10. Other toxic properties

Repeated dose toxicity. In rare cases, repeated excessive exposure to propylene glycol may cause central nervous system effects. Aspiration hazard Due to its physical properties, the substance does probably not pose any aspiration hazard.

#### 11.11. Other relevant toxicity information

Handle in accordance with good industrial hygiene and safety practice.

#### 11.12. Experience with human exposure

Health injuries are not known or expected under normal use.

# SAFETY DATA SHEET

## HEXID A4 HEAT TRANSFER FLUID

Conforming to Directive 1907/2006/EC

### SECTION 12: ECOLOGICAL INFORMATION

#### 12.1. Acute toxicity

Fish - LC50 : 40613 mg/l (Oncorhynchus mykiss; 96 h) (static test)

Daphnia and other aquatic invertebrates - LC50 : 18340 mg/l (Ceriodaphnia Dubia (water flea); 48 h) (static test)

Algae - ErC50 : 19000 mg/l (Pseudokirchneriella subcapitata (green algae); 96 h) (Growth inhibition)

Bacteria - NOEC : > 20000 mg/l (Pseudomonas putida; 18 h) Chronic toxicity

Aquatic invertebrates - NOEC : 13020 mg/l (Ceriodaphnia Dubia (water flea); 7 d) (semi-static test)

#### 12.2. Persistence and degradability

Biodegradability 81 % (anaerobic; Exposure Time: 28 d)(OECD 301 F)

Readily biodegradable 96 % (anaerobic; Exposure Time: 64 d)(OECD 306.)

#### 12.3. Bioaccumulative potential

BCF - 0.09 estimated Low bioaccumulative potential

#### 12.4. Mobility

Estimated Koc < 1, indicating very high soil mobility.

#### 12.5. PBT and vPvB assessment

Not a PBT or vPvB substance or mixture

#### 12.6. Other adverse effects

Do not flush into surface water or sanitary sewer system. Avoid subsoil penetration. This substance is not in Annex I of Regulation (EC) 2037/2000 on substances that deplete the ozone layer.

### SECTION 13: DISPOSAL CONSIDERATION

#### 13.1. Waste treatment methods

Disposal together with normal waste is not allowed. Special disposal required according to local regulations. Do not let product enter drains. Contact waste disposal services.

#### 13.2. Contaminated packaging

Empty contaminated packaging thoroughly. They can be recycled after thorough and proper cleaning. Packaging that cannot be cleaned are to be disposed of in the same manner as the product.

#### 13.3. European Waste Catalogue Number

No waste code according to the European Waste Catalogue can be assigned for this product, as the intended use dictates the assignment. The waste code is established in consultation with the regional waste disposer.

### SECTION 14: TRANSPORT INFORMATION

Not dangerous goods for ADR, RID, IMDG and IATA.

#### 14.1. EEC Regulations

**UNNO** None **Class** None **Packing Group** None

Road & Rail Transport (ADR & RID) None **IMDG** Not Applicable **ICOA** None

### SECTION 15: REGULATORY INFORMATION

**15.1** Classification Not classified as hazardous to users.

**15.2** CAS No. 57556

**15.3** Risk or Safety phrases None

**15.4** Labelling None

### SECTION 16: OTHER INFORMATION

Key literature references and sources for data taken from supplier information and data from the "Database of registered substances" of the European Chemicals Agency (ECHA) were used to create this safety data sheet. Other information - The information provided in this Safety Data Sheet is correct to our knowledge at the date of its revision. The information given only describes the products with regard to safety arrangements and is not to be considered as a warranty or quality specification and does not constitute a legal relationship.

The information contained in this Safety Data Sheet relates only to the specific material designated and may not be valid for such material used in combination with any other material or in any process, unless specified in the text.